BioHPP®

The reference for physiological substructure materials



physiological - esthetic - biocompatible

The physiological solution

Transfer of forces to the jawbone

Extremely rigid prosthetic materials withstand natural torsion of the jawbone. If they are connected to a rigid bridge structure (metal, zirconium, etc.) in the area of the pre-molars and molars, the tensile and compressive forces are increased in the area of the roots. Natural teeth are able to partly compensate for these forces but fixed osseointegrated implants do not provide this compensating effect. These forces act on implants and the bone at an unfavorable angle and, in the macro area, they also affect the physiological movement pattern and have negative effects on the dorsocranial movement capacity, CMD, osseointegration or bone atrophy.

Natural tooth


Rehabilitation with BioHPP reduces the stress caused by natural forces and the forces attributed to the prosthetic restoration.


Rigid materials prevent natural movement and transfer forces to the implants and the jawbones.

Compared to titanium, zirconium or ceramic, rehabilitation with BioHPP significantly reduces peak masticatory forces both for vertical and lateral movement. This property produces a positive effect for the patient and increases the durability of the restoration.


Transfer of peak masticatory forces to the jawbone


The Sharpey's fibers support retention of the tooth within the jaw and also absorb forces applied to the tooth. This property is no longer available after placement of an implant or for a root-treated tooth, so that peak masticatory forces are transferred directly and entirely into the jaw. From a mechanical point of view, this affects osseointegration and produces adverse effects for opposing teeth. These peak forces can be considerably reduced with abutments made of BioHPP. This is particularly important for cases of immediate restoration to ensure reliable osseointegration.

The Sharpey's fibers reduce peak masticatory forces.

To some extent, BioHPP adopts the effect of the missing Sharpey's fibers.


No damping. A restoration made of rigid materials (titanium or zirconium) provides 1:1 transfer of forces to the implant and the jawbone.

Stages of success

- ➔ Use of PEEK high-performance polymer in human medicine over more than 35 years
- ➔ More than 10 years of experience in the field of prosthetic rehabilitation
- ➔ More than 3,500 patients with BioHPP restorations
- ➔ More than 1,800 dental laboratories use BioHPP
- ➔ Documented by more than 30 university and clinical studies

See selection of studies 1-19
... on page 13

1988 the PEEK material is approved for oral use in dentistry.

2004 bredent is the first company worldwide to successfully launch a PEEK-based substructure material in the dental sector. This material is called BioXS and is still preferably used in the flasking technique.

1980

1985

1990

1995

2000

2005

1980 PEEK is an established denture resin in human orthopedic surgery due to its physiological properties and stability.

2007 bredent presents the visio.lign product line that completes the system of physiological substructure materials with esthetic materials for permanent implant prosthetics (veneers, teeth and composite) and as an alternative to ceramics.

2008 visio.link is the first and so far the only bonding agent to produce bond strength values between the visio.lign veneering products and substructure materials, such as BioHPP, that could not be achieved with any other materials and bonding agent.

2011 After completion of an international development phase, BioHPP is approved for use as universal, tooth-colored substructure material for fixed, removable, implant-supported and permanent dental restorations.

2014 The use of breCAM.BioHPP and breCAM.HIPC allows to combine the digital and analog advantages. With the 2 in 1 concept, a 24-hour overnight service can be offered to fabricate a tooth-colored substructure and a veneer with a natural appearance in a single workflow.

2010

2013 BioHPP elegance enables bredent to offer the first individual hybrid abutment (no adhesive gap) as a physiological and technically perfected alternative to titanium abutments for more than 9 implant systems now.

2015

2016 BioHPP elegance prefab is presented on the market to enable the fabrication of individual, physiological hybrid abutments (no adhesive gap) in a digital process in less than 15 minutes.

Indications

BioHPP was continuously advanced so that specific materials are available for a wide range of indications. Some examples:


Fixed restorations

- Single crowns
- Bridges (max. two pontics)
- Adhesive bridges (Maryland)

Removable restorations

- Superstructures with or without friction elements
- Secondary parts for the telescopic crown technique and bar superstructures
- Primary crowns

Implant prosthetics

- Individual abutments for 9 different implant systems (BioHPP elegance)
- Crowns and bridges (screw-retained or cemented)
- Crown and bridge substructures
- Removable superstructures
- Crown abutments
- Primary parts
- Toronto Bridge

visio.lign®

BioHPP

BioHPP
SKY
elegance

Photos: Sebastian Schuldes (MDT), Eisenach, Germany

What is BioHPP?

Biocompatible High Performance Polymer

From PEEK to BioHPP

PEEK has been used as an implant material in human medicine (finger prostheses, intervertebral discs and hip joint prostheses) for 35 years. The benefits can be attributed to the highly compatible material properties which make fusion with the bone possible. Moreover, the mechanical material properties are very similar to that of bone material.


The mechanical values of pure PEEK, however, are not sufficient for the use in more comprehensive indications and the more stringent requirements for use in the oral cavity. The base material had to be refined.

BioHPP is a partially crystalline, thermoplastic high-performance polymer based on PEEK (polyether ether ketone) filled with inorganic microparticles with a diameter of less than 0.5 µm. Physiological elasticity was retained and the combination of perfect rigidity and exceptional polishing characteristics could be achieved.

As a result, BioHPP is the only material to achieve perfect balance between:


- ✓ Elasticity and rigidity
- ✓ Weight and fracture strength
- ✓ Physiology und resistance to plaque

BioHPP is the strongest non-metallic material in dentistry

Structural formula of a PEEK molecule. The white cloud is an indicator of the ceramic filler which accounts for excellent mechanical material properties, especially for the use in dental techniques.

Competence in the field of resins

Material	BioHPP	Nature (Reference)	Pure PEEK	PMMA	PM alloy	Titanium	Zirconium
Specific weight	1.4 g/cm ³		1.3 g/cm ³	1.18 g/cm ³	19.3 g/cm ³	4.5 g/cm ³	6.5 g/cm ³
Hardness	30 HV = 294 N/mm ²		20 HV	18 HV	190 – 240 HV	300 – 400 HV	1.200 HV
Modulus of elasticity	4.200 – 4.800 MPa*	Jaw bone 2.000 – 12.000 MPa	3,600 MPa	3,600 MPa	60,000 – 130,000 MPa too hard	115,000 MPa too hard	205,000 MPa too hard
Water absorption of composites	6.5 µg/mm ³		5 µg/mm ³	19 µg/mm ³			
Water solubility	< 0.03 µg/mm ³		0.05 mg/mm ³	1-1.4 mg/mm ³	insoluble	insoluble	insoluble
Flexural strength	180 – 185 MPa		165 – 170 MPa	95 – 105 MPa			100 – 180 MPa
Bond strength (with veneering material)	> 38,8 MPa ¹⁾		20 MPa (composite)		20 – 30 MPa (with ceramic)	> 25 MPa	> 25 MPa
Thermal conductivity	low	low	low	low	high	high	low
Surface polishing characteristics	< 0,02 µm very good		poor	< 0.05 µm good	good	poor	good

* Depending on the type of processing, pressing / milling
¹⁾ When using visio.link and combo.lign opaquer

Safe

BioHPP features exclusively material properties that are identical to those of previously established materials. The use of BioHPP does not include any risks or hazards. BioHPP exhibits better values in a lot of disciplines and hence is the material of choice for the majority of indications.

Physiological

- resilient
- shock-absorbing
- non-abrasive for the residual teeth
- elastic like jaw bone
- resistant to fracture & torsion
- tolerant
- friction values are retained

Biocompatible

- metal-free
- hypoallergenic
- insoluble in water
- resistant to plaque
- prevents electrolytic conduction
- keeps inherent strength in contrast to ceramic
- no age-related degeneration
- resistant to gamma and x-ray radiation
- chemically stable

Impresses patients

- ✓ natural esthetics
- ✓ natural chewing capability
- ✓ natural feeling in the mouth
- ✓ light material
- ✓ more pleasant than rigid restorations
- ✓ optimized price-performance ratio
- ✓ does not stain

Benefits for processing

- ✓ desired shape can be achieved with limited effort (pressing, CAD/CAM, finished parts)
- ✓ easy to grind (mill), also in situ
- ✓ easy to polish
- ✓ easy to veneer

What does BioHPP offer?

Mechanical benefits

The perfect combination of elasticity and rigidity

- **Shock-absorbing** - as an abutment it supports osseointegration of implants, is suitable for immediate restorations
- **Torsionally flexible** like healthy bone, allows natural, physiological torsion of jaw braces
- **Resistant to fracture** - studies document the suitability of BioHPP for the fabrication of large-sized bridge structures with spans of up to 16 mm ^{a)b)}
- **Perfectly suitable for veneering** - BioHPP and visio.lign produce bond strength values that exceed the highest values of ceramic materials ^{c)}
- **Easy to process** - BioHPP can even be ground and polished in the oral cavity without the risk of reducing the quality of the material structure.

Prosthetic stability for permanent restorations - even for implant prosthetics

Thanks to the bound microceramic fillers, BioHPP is easier to veneer, more rigid and stable, more resistant to fracture and features better polishing properties and machinability (CAD/CAM) than pure PEEK. Moreover, BioHPP also allows the fabrication of physiological restorations ^{d)}.


Test setup with standardized samples from the University of Regensburg ^{e)}

Studies

- a) cf. study 4 University of Munich
- b) cf. study 5 University of Munich
- c) cf. study 8 University of Regensburg
- d) cf. study 2 University of Jena
- e) cf. study 17 University of Regensburg

... on page 13

Transfer of peak masticatory forces to the jaw bone

BioHPP compensates the lack of natural shock absorption by the Sharpey's fibers. The elasticity of BioHPP reduces the peak masticatory forces and directs the forces to the bone via the implant over an extended period. As a result, immediate rehabilitation with implants and optimized osseointegration are enabled.

Physiological elasticity

The modulus of elasticity of BioHPP that corresponds to the average value of the jaw bone differs by a factor of 27 from the value of rigid substructure materials, such as titanium and ZrO₂, which are still frequently used.

E-modulus	Factor			
Jawbone	BioHPP	Gold	Titanium	Zirconium
1,000 - 12,000 = 4,200 - 4,800 MPa	=	x 20	x 25	x 27

Fault-tolerant, versatile processing

Restorations made of BioHPP - large-span bridges or single, individual abutments - can be fabricated in the thermo-plastic pressing technique or in a CAM process. Reworking using grinding and polishing processes can be easily and quickly completed and does not compromise the material quality unlike for ZrO₂. It can and should be carried out in the oral cavity.


Biological benefits

Superior physiology from the biological view

- **Biocompatible** – BioHPP is a class 2a medical device and complies with all relevant DIN standards, is not cytotoxic and hence biocompatible (DIN 10993-05,10, 11, 03, 12)
- **Metal-free**, no ion exchange, not allergies, no metal taste ^{d)}
- **Plaque-neutral 1**, the official polishing protocol for BioHPP shows a surface smoothness which partly exceeds the one of a successful veneer^{d)} This can be attributed to the small particle size of $< 5\mu\text{m}$ of the bound inorganic microparticles
- **Plaque-neutral 2**, the low water absorption of $6.5\ \mu\text{g}/\text{mm}^3$ avoids plaque accumulation and hence odor generation and discoloration
- **Kind to the gingiva** – soft tissue management is considerably facilitated by the high level of acceptance and integration of BioHPP in particular since individual BioHPP abutments are suitable one-time abutments in immediate restorations^{g)}

Kind to the gingiva


Apparent after removal of the gingiva former: a bleeding margin indicates significantly improved gingival attachment compared to the titanium abutment.

Photos from the Department for Implantology of the Dental Clinic IRCCS Fondazione Policlinico Ospedale Maggiore Università degli Studi di Milano (Scientific & Clinical Cases: „Physiologische Prothetik“, see page 16).

Ease of polishing and cleaning

BioHPP can also be easily cleaned by the patient using a soft toothbrush without roughening the surface. Professional and chairside tooth cleaning can be conveniently performed with a soft abrasive agent and subsequent polishing is carried out with standard tools. Surface roughness of $0.05\ \mu\text{m}$ is achieved to protect against discoloration and plaque accumulation. In the direct comparison, BioHPP provides better polishing properties than successful dentures and veneering composites.^{h)}

Homogeneous and fine-grained surface of BioHPP at a magnification of 1000x under the electron microscope.


Studies

- f) cf. study 17 University Clinic of Tübingen
- g) cf. Scientific & Clinical Cases: „Physiologische Prothetik“
- h) cf. study 7 University Clinic of Cologne
... on page 13

Immediate restorations – clinical advantages of BioHPP (SKY) elegance

Conclusion

The analysis of the clinical cases and animal experiments document that – in cases of immediate loading – the new individual abutment BioHPP SKY elegance is a reliable alternative to industrial titanium abutments. Moreover improved healing of the soft tissue is achieved.

Histological analysis of the SKY elegance abutments. Details on the platform switch and the gingival attachment:

- a) After 4 weeks
 - b) After 8 weeks
- Prof. Dr. José Eduardo Maté Sanchez de Val, PhD, MSc, DDS, Murcia, Spain; this study will be published .

What does BioHPP offer?

Prosthetic benefits

Outstanding quality of the restoration

- Stability - thanks to the very high mechanical strength values, BioHPP is suitable for very large bridge substructures (spans of up to 16 mm) and for removable restorations in implant prosthetics^{a)}
- Resistant to abrasion - telescopic restorations made of BioHPP demonstrate exceptional durability of the friction function^{b)}
- Can be perfectly veneered - BioHPP in combination with the visio.lign system produces higher bond strength values than restorations made of NPM (veneered with ceramics) and ZrO₂^{a)}
- Light - patient acceptance is increased thanks to the light weight of the restoration; BioHPP substructures are 4 times lighter than ZrO₂ substructures
- Tooth-colored - exposed areas made of BioHPP (available in white and tooth-colored) can hardly be noticed, dark areas do not need to be compensated by the veneer
- Natural feeling in the mouth, weight, thermal conductivity, elasticity, smoothness and the physiological integration into the body system ensure that the patient forgets that he is wearing a denture

Studies

a) cf. study x University of Munich
 b) cf. study 12 University Clinic of Cologne
 c) cf. study 2 University of Jena

... on page 13


Wearing comfort

Patients who are able to make a comparison prefer restorations made of BioHPP (instead of ZrO₂, for example) because of the light weight of the denture and the natural feeling in the mouth.

Weight	Factor			
Jaw bone	BioHPP	Gold	Titanium	Zirconium
1.3 - 1.4 g/cm ³	x 1	x 14	x 3	x 5

Reliable friction


Compared to ZrO₂ and non-precious metal, BioHPP reveals superior consistency in wear tests. Based on the resistance to abrasion, the University Clinic of Cologne even recommends: „The use of BioHPP primary elements with BioHPP secondary elements should be preferred.“^{b)}

Result diagram of the University Clinic of Cologne^{b)}

Reliable bonding to the veneer

Studies have shown that BioHPP in combination with the visio.lign system produces the highest bond strength values compared to silanized CoCr or metal veneered with ceramic. The advantage of the visio.lign system lies in the fact that it contributes to the physiological benefits of BioHPP. This combination ensures perfect, most efficient and economic restorations.

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and numerous others; please request information!

How is BioHPP cleaned and polished?

By the patient

The patients should preferably clean his teeth every day using a soft to medium hard toothbrush. The use of rotary electric toothbrushes is recommended; ultrasonic toothbrushes, however, should not be used since they may have adverse effects on the surface.

In the dental practice

Professional tooth cleaning

Restorations made of BioHPP can be easily cleaned and subsequently polished in the dental practice using standard instruments and materials for soft high performance polymers. For more and detailed information, the „BioHPP Praxis Guide“ is available.

Sterilization

The individual and adhesive gap-free elegance abutments can be sterilized using steam sterilization (autoclave) in the vacuum procedure. For this purpose 3-times fractionated prevacuum must be generated with a sterilization time of 4 minutes and a temperature of 134°C +/- 1 °C.

Dental hygienist Ms Vesna Braun

How is BioHPP bonded?

Type of bonding	BioHPP® crowns and bridges to...	Bonding systems				
		Metal/alloy abutments	Zirconium dioxide abutments	BioHPP® abutments	Hard tooth substance (dentine/enamel)	Use visio.link on BioHPP®
permanent	Adhesive – with conditioning/primer using composite cement, e.g. Panavia F 2.0 (Kuraray), Variolink II (Ivoclar), NX-3 (Kerr)	✓	✓	✓	✓	✓
	Self-adhesive composite cement 110 µm (sand-blasting), e.g. Rely X Unicem (3M Espe)	✓	✓	✓	●	●
	Glass ionomer cement, e.g. Ketac Cem (3M Espe)	●*	●*	●	●*	X
	Zinc phosphate cement (e.g. Harvard)	●	●	●	●*	X
temporary	Zinc oxide, eugenol-free cement (Tempbond by Kerr)	✓	✓	✓	●*	X
	A-Silicone-based fixing cement (Temposil 2 by Coltène Whaledent)	✓	✓	✓	✓	X

* Use only for preparation angles up to 5°

Type of bonding	BioHPP abutment with framework materials made from...	Bonding systems				
		use visio.link on BioHPP®	Dental alloys	Zirconium dioxide	BioHPP®	e.max (lithium disilicate/lithium silicate) silanized
permanent	Adhesive – with conditioning/primer using composite fixing cement, e.g. Panavia F 2.0 (Kuraray), Variolink II (Ivoclar), NX-3 (Kerr)	✓	✓	✓	✓	K
	Self-adhesive composite cement 110 µm (sand-blasting), e.g. Rely X Unicem (3M Espe)	●	✓	✓	✓	X
	Glass ionomer cement, e.g. Ketac Cem (3M Espe)	X	●*	●*	●	X
	Zinc phosphate cement (e.g. Harvard)	X	●	●	●	X
temporary	Zinc oxide, eugenol-free cement (Tempbond by Kerr)	X	✓*	✓*	●	X
	A-Silicone-based temporary cement (Temposil 2 by Coltène Whaledent)	X	✓	✓	✓	X


* Use only for preparation angles up to 5°

✓ = optimum K = only to be used for crowns ● = generally possible X = not recommended

Scientific & Clinical Cases

Physiological Prosthetic Immediate single-tooth restoration

Various practical cases, both scientifically and clinically-documented and complete with illustrations. Discover new approaches to restoration options and obtain ideas you can use in your own laboratory.

Available in German REF 9929760D
and English REF 992976GB


Available in German REF 9929770D
and English REF 992977GB

Scientific & Clinical Cases online

The online version of Scientific & Clinical Cases can be viewed by scanning the QR code or by visiting www.bredent-medical.com/en/scientific

Other offers that may be of interest to you:

REF 000535GB

REF 000722EX

REF 000588GB

000547GB-20160324 Mistake and subject to change reserved

DENTAL INNOVATIONS
SINCE 1974

breident
group

breident GmbH & Co. KG
breident medical GmbH & Co. KG

Weissenhorner Str. 2 · 89250 Senden · Germany
T: +49 7309 872-441 · www.bredent.com · @: info@bredent.com