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FAST  
PROTEC

**XPdent** Corp.

Sales | 1.877.328.3965

Tech Support | 1.877.329.3965

Web | [www.xpdent.com](http://www.xpdent.com)

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FAST  
PROTEC

# THE ESS

## PART OF YOUR

# ESSENTIAL DAILY WORK

“

 Fast Protec has an essential nature, since it synthesizes long and expensive traditional processing in a fast, clean and easy way of working. We invest in technology and regular research for constant innovation as aesthetics, ergonomics and functionality play decisive roles in reaching simplicity. As our approach always starts with dental technicians needs, we created a complete system, made up of flasks, articulators, verticulators, resins and silicones to perform all dental processing. This is Fast Protec: precise protocols that will guide you step by step to make your workflow faster and more organized.

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# Flasks Muflas

## ALU BIG

🇺🇸 Alu Big is a flask made of aluminium with an ergonomic design. By using silicone templates, it is possible to pour, inject or press the resin to obtain, in a fast and precise way, the most part of dental processing as removable dentures, relines, removable partial dentures, crown and bridges, surgical templates and implant dentures.

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---

## IMPLANT BOX

🇺🇸 Implant Box is a flask with a minimal and compact design that allows you to cast the model in plaster directly on the base. Thanks to the integrated split cast it is possible to perfectly fit the base on Fast Protec Articulator, Verarte and All-In. It is provided with two tops: the metallic one is useful to perform crown and bridges and implant dentures whereas the clear one is perfect to work with light-curing composites.

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## BIG AND SMALL GLASS

🇺🇸 Big and Small Glass are transparent flasks with a POM base and two clear lids held by secure locking screws to ensure optimal results and no increases in the vertical dimension. The sizes of these flasks make them compatible with any polymerization chamber available on the market, and the double lids have been designed for the duplication of crown and bridge wax-up in light-curing composite.

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## LIGHT GLASS

🇺🇸 Wider than traditional transparent flasks and more stable with its three screws, Light Glass allows dental technicians to perform full implant dentures working simultaneously on composite teeth and gingiva thanks to the provided three lids and the removable three pins. Its completely clear structure guarantees a fast and total light-curing of the elements obtaining perfect and precise details.

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# Articulator and Verticulators

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### ARTICULATOR

 The unique Articulator of this system is provided with magnetic condyles fixed at the standard 30 degrees to perfectly reproduce all mouth movements. Diagnosis and treatment planning are simplified thanks to the posterior incisal guide pin that gives total visibility at 360°. Functional and easy to use, Articulator enables to set the model directly in articulation through the use of silicone instead of plaster.

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---

### VERARTE

 Verticator and Articulator combined in one unit in a revolutionary design. The posterior rod pin guarantees a perfect and complete visibility of the prosthesis and simply by changing the upper support of the Articulator, it is possible to transform it in a Verticator to perform ceramic moulding or plan implant dentures.

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## ALL IN

 All-In is a combination of Verticulator and Flask composed by one base with three rod pins and two interchangeable lids. The equilateral base enables the rotation of the model and the template at 360°, allowing the dental technician to have the maximum visibility of his work. Essential item for the pre-visualization of the case and the casework planning as a Verticulator, it could be used for the injection or pouring technique as a Flask.

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## Implatologia

### EASY BAR

 Easy Bar is a kit of stainless steel preformed bars that facilitate the realization of a stable, precise and reinforced structure for all kind of fixed dental prosthetics in just 6 hours. Thanks to their different sizes (from 6 to 30 mm), these bars could be adapted by bending to any position of the implants and are compatible to every implant abutment on the market.

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---

### JIG PLATFORM

 Jig Platform is ideal for positioning and aligning drill holes for abutment holes in dentures. It is provided with a secure magnetic attachment on the base and an integrated split cast.

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# Accessories Accesorios

## DOUBLE BASE

🇺🇸 Two base formers made of aluminium with rubber cuffs useful to cast the model in plaster. Thanks to the split cast design on the base, the model becomes compatible with any product of the Fast Protec System, fitting in a precise and clean way.

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---

## ALU BIG CLEAR TOP

🇺🇸 This Clear Top transforms Alu Big into a transparent flask to work with light-curing composites.

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## TEETH CLEAN

 Basket with 14 individual compartments for steam cleaning teeth securing each denture tooth position.

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---

## MODEL PLATE

 Model Plate is a positioning device for setting edentulous models on the articulator.

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## CLEAN TEMP

 Split cast base for temporaries crowns and bridges. Minimum waste of material thanks to its three undercuts.

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## BASE FOR PLASTER

 Base for adapting any articulator on the market to the Fast Protec System.

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---

## UNIVERAL CLEAN

 Universal Clean is a small base that helps fit the plaster model on the Articulator through the use of Universal Blue Silicone.

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---

## UNIVERAL BASE SPLIT CAST

 Small split cast base to adapt other bases available on market to the Fast Protec system.

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# Silicones Siliconas

## UNIVERSAL BLUE SILICONE

 Kneadable silicone for the first templates. It guarantees an extremely precise detail reproduction without distortion even at high temperatures and an excellent long-term dimensional stability.

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---

## PUTTY SILICONE

 With ideal handling characteristics and times, Putty Silicone is an economic rigid bulk filler silicone.

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50  
SHORE

## TRASPA SILICONE

 Crystal-clear silicone perfect for processing with light-curing composites.

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Content: 1 kg A + 1 kg B

## REGULAR PRECISION

 Auto-curing and precise silicone useful to reproduce temporaries caps or for temporaries restoration.

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
clear

## ACRYLIC FOR DENTURES MICRO PEARL

 Cold-curing resin for the injection, pouring or pressing techniques. Different colours available but the same stability and no shrinkage to meet all your needs: full and partial dentures, repairs and relining could be performed in an easy and fast way. Top aesthetic prosthesis with less residual monomer.

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

 Full Denture  
 Prótesis Móvil

 Immediate Loading  
 Carga Inmediata


 Temporary Crowns  
and Bridges  
 Provisionales

 Ceramic Moulding  
 Estampe de la ceramica

 Orthodontic Splint  
 Férula

 Implant Denture  
with Composite  
 Toronto en composite

# Full Denture Prótesis Móvil


Performed by: Dt. Vincenzo De Feo


## PRODUCTS USED:

- Alu Big
- Articulator
- Teeth Clean

## MATERIALS USED:

- Acrylic for Dentures Micropearl
- Universal Blue Silicone 70 Shore A+B
- Putty Silicone 70 Shore A+B
- Fast Alginate
- Gengicolors Kit

1. Wax-up
2. Occlusion check on the Articulator
3. First silicone template on the Alu Big
4. Second template
5. Teeth degreasing
6. Repositioning of the teeth on the silicone template
7. Isolation of the model
8. Closing of the flask
9. Resin pouring
10. Opening of the flask
11. Polishing
12. Final occlusion check

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
Immediate Loading processing shown in this page has been performed by Dt. Umberto Scilipoti

1. Implant denture model
2. Control of the bar position (previously fixed in mouth)
3. Model set on Articulator
4. Easy Bar structure welded
5. Easy Bar structure cemented on the abutments
6. Wax-up
7. Silicone templates
8. Opacified structure
9. Teeth insertion
10. Injection of the resin
11. Opening of the muffle
12. Final result after polishing

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# Temporary Crowns and Bridges Provisionales


Performed by: Dt. Francesco Raffo

## PRODUCTS USED:


- Alu Big
- Double Base
- Base for Temp
- Articulator

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

### MATERIALS USED:

- Universal Blue Silicone 70 Shore A+B

1. Occlusion check on the Articulator
2. Silicone template
3. Check of the spaces
4. Ceramic moulding
5. Dentine layering
6. First firing
7. Enamel layering
8. Final work after the second firing

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25min

~~60min~~

1. Occlusion check
2. Bite check on the Articulator
3. Bite modelling
4. First silicone template
5. Second silicone template
6. Flask opening
7. Isolation of the model
8. Preparation of the resin
9. Injection of the resin
10. Flask opening
11. Occlusion check on the Articulator
12. Final result after polishing

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Performed by: Dt. Luigi De Stefano

## PRODUCTS USED:


- Articulator
- Light Glass
- New Light Box

## MATERIALS USED:

- Universal Blue Silicone 70 Shore A+B
- Traspasilicone 22 shore A+B
- Traspasilicone 50 shore A+B
- Comp Metal Primer
- Comp Metal Connector
- Comp Bonding
- Comp Flow
- Intensive Comp

3h

~~10h~~

1. Framework treatment and wax reproduction of the dental anatomy.
2. Setting of the wax-up on the BASE (B) with silicone and creation of injection canals. Traspa Silicone 22 shore pouring on COVER n°1 to obtain the template of the complete wax-up.
3. Positioning of the wax-up on COVER n°1, creation of injection canals and Traspa Silicone 50 shore pouring.
4. Wax-up set on the BASE (B) and cut back.
5. Creation of injection canals for the enamel and closure of the flask with COVER n°2. Traspa Silicone 22 shore pouring.
6. Metal framework cleaned removing wax. Preparation of the metal framework using Metal Primer, Metal Connect and Comp Flow Opaquer.
7. Injection of Comp Flow Dentine on COVER n°2 and light-curing in the New Light Box.
8. Painting with Intensive Comp colours and light-curing in the New Light Box. Injection of Comp Flow Incisal enamel on COVER n°1 and light-curing in the New Light Box.
9. Injection of Gengiva Flow composite on the BASE (B) and light-curing in the New Light Box.
10. Final light-curing in the New Light Box. Finishing and polishing.

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# OUR VALUES

1

 We believe essentiality is the result of experience, research, collaboration and technology. It is a way of working with efficiency and efficacy.

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2

 Thanks to the high know-how of our team, we are able to design and produce innovative products that could facilitate your daily work. Everyone contributes with his own talent, diligence and expertise always considering the dental technician as the centre of each choice.

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3

 We listen to our customers needs to find together the best solution.

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4

 In the last 25 years the best materials, and constant innovation have always been our prerogatives to guarantee very high standards in every process.

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## QUALITY AND CUSTOMER SERVICE

 Customer is at the core of all our choices: we grant high quality products and a constant service to support you in your daily operations avoiding wasting time and material. Our customer service:

- Pre-sale and post-sale assistance
- Constant innovation of our products, materials and accessories to improve our system
- Training courses and follow-up meetings with dental experts

Quality is made up of little details. That is why we use nickel-free anodised aluminium to assure you the unlimited life of our system.

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# Polymerization units Polimerizadoras

## VICTORIA

 Compact, easy and functional polymerizing machine. It is possible to cook cold and heat-curing resins with water input and to harden silicone templates with manual air loading through button up to 6 bar pressure. Victoria is provided with heat strips that heat water gradually until a maximum temperature of 120 degrees and an automatic counter that starts decreasing time once the pre-set temperature is reached maintaining it constant.

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## NEW LIGHT BOX

 The New Light Box is a compact and modern polymerization unit provided with last generation multi-die LEDs for light-curing composite resins. A low consumption machine that works on a vacuum pump to avoid the dispersion layer, the New Light Box is perfect to light-cure the Fast Protec wide range of composites at a wave frequency between 360 and 520 nm.

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## GENGICOLORS KIT

 Specifically studied for Fast Protec acrylic micro pearled resins, Gengicolors is a kit composed by two pre-mixed light and dark cold-curing resins and five stain colours for the characterization of marginal and attached gingivas. Optimal to give dentures a natural and individualized look.

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Content: light pink powder 100 gr; dark pink powder 100 gr; liquid 100 ml; stains 10 gr


## TEMP POWDER AND LIQUID

 High-tech acrylic for long lasting crowns and bridges. Temp Powder is a cold-curing resin with a high precision and abrasion resistance that avoids plaque deposits.

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Content: dentine 100 gr; incisal 100 gr; liquid 100 ml

# Isolating Agents Aislantes

## FAST ALGINATE

 Separating agent, based on alginate, useful to isolate the plaster model before injecting or pouring the resin.

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Content: 500 ml

## SILICONE SPRAY

 Separating spray excellent to insulate silicone templates.

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Content: 400 ml

## GENGIVA FLOW

 Light-curing gingiva composite ideal for gingiva restorations and reproductions, implant suprastructures, crowns and bridges. Thanks to Fast Protec wide range of colours, it is possible to characterize full or partial dentures with incredible final aesthetic results.

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Content: 1 syringe 3g

## MASTER POLYURETHANE

 Polyurethane resin ideal for fabricating precise and resistant models as well as crowns and bridges. Flowable and easy to trim, Master Polyurethane has a dimensional stability that guarantees precise details of all works with a shrinkage percentage close to zero.

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## UNIPATT RESIN

 Self-curing castable resin for telescopes, super- and supra-structures using the brush technique.

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We would like to thank Vincenzo Cuiuli,  
inventor of the [Dental Processing System](#) for his unique contribution.

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**XPdent** Corp.

Sales | 1.877.328.3965

Tech Support | 1.877.329.3965

Web | [www.xpdent.com](http://www.xpdent.com)