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Opaqer combo.lign

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Instructions for use

Opaque combo.lign

Before using the product, please read these instructions for use carefully.

1. Indication range

1.1 Product description

Opaque combo.lign is an opaque resin (composite) for coating porcelain veneer crowns, porcelain veneer bridges and metal frameworks with the respective tooth shade of the veneer. Opaque combo.lign is a dual-harden polymer resin which is polymerized by heat (self-curing composite) and using a light source with a suitable wavelength (370 to 500 nm). 4 light-curing shade pastes in syringes of 4 g each (approx. 3 ml) and a catalyst paste in a separate syringe (cont. 4 g) are available. Opaque combo.lign is intended to be used for fixed and removable dental restorations and for extraoral processing. Opaque combo.lign may only be used in combination with an adhesive (bonding) system. All inorganic (silicization with ensuing silanization) and chemical primers to be used for metal-resin (composite) bonding are suitable.

1.2 Indication

Coating (covering) metal frameworks based on precious metal-free (PMF/NPM) and precious metal-containing (PM) alloys using a primer (adhesive or silane) for

- Telescopic and conical crowns (cover denture)
- Veneer crowns, veneer bridges
- CoCr structures
- Metal frameworks or tertiary structures

Moreover Opaque combo.lign is also suitable for coloring metal-free zirconium dioxide frameworks and thermoplastic framework material BioXS which are subsequently veneered using novolign veneers or combo.lign and crea.lign veneering materials.

1.3 Contraindication

Opaque combo.lign is not suitable for intraoral processing.

1.4 Shades and classifications for the classic A-D shade system

Light:	A1 / A2 / A3 / B2
Medium:	A3,5 / B3 / D3
Intensiv:	A4 / C3
GUM:	Gingiva

2. Preparatory steps

2.1 Conditioning of metal frameworks

The metal framework is sandblasted with 110 µm aluminium oxide at a pressure of 3 to 4 bars (surface enlargement). Then the surface is wetted with Primer II from GC (bonding agent for bonding composite to metal) and cleaned with steam (Note: adhere to the processing instructions for Metall Primer II!) until the surface is dry. Alternatively, silanization can be carried out using inorganic procedures (for example Silano-Pen or Rocatec); please adhere to the instructions for use!

2.2 Conditioning of zirconium dioxide frameworks and crowns

The framework or crown is sandblasted with aluminium oxide (50 to 110 µm) at a max. pressure of 2 bar (activation). Then Silano-Pen (REF 320 0047 0) is used for silicization and silanization (bonding agent) is carried out in accordance with the instructions for use for Silano-Pen (REF 320 0047 0).

2.3 Conditioning of BioXS

Sandblast the composite framework or crown 110 µm aluminium oxide at a pressure of 2 to 3 bars (to enlarge the surface). Then wet with light-curing PMMA & Composite Primer visio.link (REF VLP/PMMA10) and polymerize subsequently in the light-curing unit (for example 90 seconds in the Heraeus Kulzer UniXS unit) in accordance with the processing instructions for visio.link.

3. Processing

3.1 Mixing, applying and polymerizing

After conditioning the framework, add Opaque combo.lign to one part of light-curing shade paste and one part of catalyst paste on a mixing plate, tray or block and stir thoroughly using a plastic spatula. Use a brush to apply the mixed paste to the framework; processing time span is 5 to 6 min. at a room temperature of 20° C (extending the processing time, see 3.2). To achieve optimal bonding strength, hardening with a light-curing unit is required, for example Uni XS 2 x 180 sec. Do not use polymerized (by chemical hardening) paste or residues left on the mixing plate, tray or block.

Information on formation of inhibition layer!
Thanks to dual-hardening, a slight larger inhibition layer (sticky layer) is formed, which is required for ensuring bonding of additional layers or luting of the veneers. Therefore Opaque combo.lign must not be hardened or polymerized under vacuum. The inhibition layer does not indicate insufficient hardening; the layer must not be removed!!

3.2 Adjusting the processing time:

The processing time span can be adjusted to approx. 5 to 6 min for processing temperatures of e.g. 25° C by changing the mixing ratio with a maximum of two parts of light-curing shade paste and one part of catalyst paste (2:1). If the ambient temperature is increased by 10° C, the reaction rate will double!

3.3 Light polymerization units and polymerization times

Different light sources are used in light polymerization units (dental laboratories and practices). Light sources are required which cover a wavelength range of 370 to 500 nm. LED light polymerization units with wavelengths of 320 to 400 nm are not suitable. Suitable polymerization units and polymerization times are available at request.

4. Durability and storage

The storage in the comprehensive system can be extended by separating the light-curing shade paste and the catalyst paste. If proper storage of Opaque combo.lign at 24° C is ensured, the light-curing paste of Opaque combo.lign features a shelf life of 3 years. The catalyst paste contains benzoyl peroxide which may cause allergic reactions. Sensitization may be caused in susceptible persons. In case of allergic reactions or known allergies to any of the constituents, such as resins based on (meth-) acrylates, the product should not be used. Avoid skin contact! In case of contact with eyes or extended contact with the oral mucosa, rinse immediately with copious amounts of water.

5. Safety recommendations and hazard warning

Opaque combo.lign is composed of acrylate oligomers, fillers (for example dental glass powder), catalysts and color pigments. The catalyst paste contains benzoyl peroxide which may cause allergic reactions. Sensitization may be caused in susceptible persons. In case of allergic reactions or known allergies to any of the constituents, such as resins based on (meth-) acrylates, the product should not be used. Avoid skin contact! In case of contact with eyes or extended contact with the oral mucosa, rinse immediately with copious amounts of water.

6. Additional information

The information contained in these instructions for use is always updated according to the latest knowledge and experience. Therefore we recommend reading the instructions for use again before using a new package. These non-binding user recommendations are based on our own experience. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since brendent does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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