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crea.lign

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Instructions for use

crea.lign

Before using the product, please read these instructions for use carefully.

1. Indication range

crea.lign is a light-curing composite for completion of composite veneers on metal and zirconium dioxide frameworks. crea.lign is especially used for customizing, building up and finishing veneers with prefabricated veneers (novo.lign A and novo.lign P) and teeth made of PMMA or composite materials. crea.lign is suitable for oral and extraoral processing of fixed and removable restorations and hence also suitable for repairs.

1.1 Product description

crea.lign features a fluid consistency, is filled exclusively with inorganic microparticles and is therefore elastic and can be easily polished. crea.lign is polymerized using polymerization units or dental hand-held curing devices at a wavelength of 370 to 500 nm. crea.lign is supplied in syringes containing 4 g each. To ensure the maximum processing time span of the light-curing material, it is recommended to exclusively use black plastic adapters (REF 320 0094 0) which avoid the penetration of light.

1.2 Indication

- a) Completion (finishing) or repair of veneers with prefabricated veneers or teeth made of PMMA or composite materials on
- Telescopic and conical crowns (cover denture)
- CoCr structures
- Veneer crowns and bridges made of metal and zirconium dioxide
- Implant-supported structures
- Temporary restorations
- b) Building up the palatal and lingual areas for completion of veneers
- c) Customizing or repairing veneers

1.3 Contraindication

Occlusal contouring in the posterior region, customizing of veneers and teeth made of PMMA on vestibular surfaces, bruxism.

1.4 Shades

crea.lign is available for incisors, dentine and gingiva. The shade and designation of dentine shades is based on the classic A-D shades and matched with the novo.lign A and novo.lign P veneers and combo.lign.

2. Processing

2.1 Conditioning of metal frameworks

The metal framework is sandblasted with 110 µm aluminium oxide at a pressure of 3 to 4 bars (surface enlargement). Then the surface is wetted with Primer II from GC (bonding agent for bonding metal / opaque); allow the primer to evaporate subsequently (adhere to the processing instructions for the bonding agent!). Alternatively, silanizing using inorganic procedures, e.g. Silano Pen or Rocatec, can be carried out. After conditioning the metal, apply the respective shade of opaque material, e.g. Gradia opaque from GC (light-cured composite system) and polymerize (adhere to the processing instructions for the opaque and the light-cured composite!).

2.2 Conditioning of zirconium dioxide frameworks

Sway the Silano Pen across the entire surface to be veneered for approx. 5 sec. Once the flame-treated surface has cooled down, apply bonding agent and let it dry for 3 min. Process and polymerize the dentine material in accordance with the manufacturer's instructions.

2.3 Conditioning the veneer / teeth or composite veneer

Roughen the veneer or the base of the composite veneer by sandblasting with aluminium oxide (110 µm) at 2 to 3 bar or using tungsten carbide burs (surface enlargement). Then wet with light-curing PMMA & composite Primer visio.link, (REF VLP MMA10), and polymerize using a light-curing unit (e.g. 90 seconds in UniXS, Hereaus Kulzer) in accordance with the instructions for use for visio.link.

2.4 Application of crea.lign

The selected shade of crea.lign is applied to the conditioned veneer or base (combo.lign). A maximum layer thickness of 1 mm for each layering or polymerizing process must be ensured. Then polymerization is carried out using a light-curing unit. It is recommended to polymerize for 360 seconds (2 x 180 seconds) using light-curing units such as DencatolcolorXS or UniXS from Hereaus Kulzer (information about other polymerization units, see section 2.5). If several layers are applied, intermediate polymerization of at least 180 seconds (UniXS) is required.

Processing information:  
Composites such as combo.lign and crea.lign form an inhibition layer during polymerization. After grinding (milling) combo.lign, the ground surface must be conditioned with visio.link (conditioning is carried out as described under item 2.3) before additional material is applied. After the application of the final layer of crea.lign, polymerization must be carried out immediately. To avoid inhibition layers, light-curing units featuring equipment/accessories to reduce or avoid the formation of oxygen or gels that reduce the inhibition layer may be used. Prior to polishing, the gel or the inhibition layer must be removed with steam or isopropanol.

2.5 Light-curing units and polymerization times

Depending on the application, different light-curing units are used in dental laboratories and practices. To achieve complete curing (hardening) of crea.lign, a wavelength range of 370 to 500 nm is required. Pure UV units with wavelengths up to 400 nm are not suitable. Your bredent sales representative will recommend suitable units.

2.6 Surface treatment and polishing

Cross-cut tungsten carbide burs with different shapes are suitable for optimal surface treatment. Interproximal spaces are finished using fine separating discs such as Transflex, (REF 340 0005 0). Pumice polishing paste, (REF 520 0016 0) and a soft goat-hair brush are recommended for profiled surfaces that are to be smoothed slightly. Further polishing is done using soft goat-hair brushes and the Acrypol polishing pastes, (REF 520 0017 0). Use a leather buff, (REF 350 0066 0), and Abraso-Starglanz, (REF 520 0016 3), for high luster polishing.

3. Durability and storage

crea.lign must be stored at 4 - 24° C and protected against heat and sunlight. crea.lign must not be used after the expiry date.

4. Safety recommendations and hazard warning

crea.lign is composed of various acrylate oligomers, fillers (e.g. silanized dental inorganic fillers), catalysts and color pigments. Sensitization may be caused in susceptible persons. In case of allergic reactions or known allergies to any of the constituents, such as resins based on (meth-) acrylates, the product should not be used. Avoid skin contact! In case of contact with eyes or extended contact with the oral mucosa, rinse immediately with copious amount of water and consult a doctor.

5. Additional information

The information contained in these instructions for use is always updated according to the latest knowledge and experience. Therefore we recommend to read the instructions for use again before using a new package. These non-binding user recommendations are based on our own experience. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since bredent does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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