

# crea.lign®

Veneering composite



*Implant-prosthetic restoration*

© 2016 by DT Andreas Lüdtke,  
Bayreuth, Germany



## create natural beauty

Red-White Aesthetics

powered by  
**visio.lign®**



## visio.lign®

The aesthetic and functional system with unique material options, helping you to create natural beauty, function and physiology. The open system provides flexibility and freedom thanks to six components designed to compliment one another; crea.lign®, novo.lign®, neo.lign®, visio.CAM, visio.paint & bond.lign.

The combination of these components allows a variety of workflows and implant prosthetic solutions, each with perfect results. visio.lign® - the system for stability and reliability. The result: durable prosthetic restorations featuring colour stability and superior load capacity.


Truly natural and durable

**crea.lign®**

## Δ E - tendency to discolouration

Colouring substances:

Storage in coffee, tea, tobacco, red wine




New formula for crea.lign®

Gel and GUM (light, pink, rosa, red, lila).

Benefits:

- Durable restorations thanks to high fracture toughness
- No fissures or cracks thanks to elasticity on the level of denture resins
- Unsurpassed curing depth thanks to special light transmission
- Direct layering on BioHPP® is possible, no opaque required (when using visio.link®)
- Natural effect without translucency in the GUM materials
- Convenient modelling thanks to the improved stability of crea.lign® gel
- Simple application of the pastes thanks to the enhanced consistency and reduced tackiness

## “Unlimited flexibility”



Red & White with crea.lign®.

Fast. Easy. Durable. Relinable. Extendable.

The solution for:

- implant-supported, fixed or removable permanent restorations,
- natural customisation of full dentures,
- all substructure materials such as BioHPP® and PEEK,
- complete coating even of full-arch restorations (e. g. Prettau bridges)

**"Your selection ...**



Layering with  
**crea.lign® Gel**

**... of technique"**



Layering with

**crea.lign® Paste**

## Application example

8



crea.lign® Paste PO (Paste Opaque) is applied in the mucolabial fold and crea.lign® Paste PL (Paste Light) is applied to the gingival margin. Then the surface is "stippled".

crea.lign® red/lila is filled in the resulting recesses and polymerised.

Finally the surface is sealed with crea.lign® Transpa, smoothed and polymerised (end polymerisation).



Photos:

© 2016 by DT Andreas Lüdtke, Bayreuth, Germany



## crea.lign® Paste


“ The outstanding modelling properties and natural shade are what I appreciate about crea.lign® GUM Paste. Work can be easily completed and the results are fascinating. ”



**DT Andreas Lüdtke**

## crea.lign® Gel





“ The crea.lign® Gel GUM veneering material takes the benefits provided by substructures made of high-performance polymers to the next level. Feeling in the mouth, aesthetic appearance and durability are a unique combination. ”



**MDT Philipp von der Osten**

## Red-White Aesthetic Set

Individualise your restorations to achieve high-end aesthetics using the crea.lign® GUM materials in the visio.lign® system. Offer your customers this additional individual service and increase the profitability in

your laboratory. Use your dentures to increase the reputation of your laboratory.

REF CLIGNSETG



### Contents:


- |  |  |
|--|--|
| 5 x 5 g crea.lign® Gel, GUM (light, pink, rosa, red, lila) | 1 x Unique Brush, size 4                             |
| 2 x 3 g crea.lign® Paste, GUM (PO, PL)                     | 1 x One-hand use handle with UV LED (crea.lign® Gel) |
| 1 x 5 g crea.lign® Gel, Modifier beige                     | 1 x Mixing block                                     |
| 1 x crea.lign® Gel, Transpa clear                          | 1 x Pin Cover + application cannulas                 |
| 1 x 4 g crea.lign® Opaque, GUM                             | 2 x Rotary piston (crea.lign® Paste)                 |
| 1 x 10 ml crea.lign® Modelling Liquid                      | 1 x Layering instructions                            |
| 1 x 10 ml visio.link®                                      |  |

## crea.lign® Gel or Paste

### Dentin



Gel



Paste



## crea.lign® Gel

### GUM



Gel



## crea.lign® Paste

### GUM



Paste



## crea.lign® Gel

### Enamel & Incisal



Gel



### Modifier



Gel



### Transpa Clear



Gel



# visio.lign®

## The Aesthetic and Functional system

### Overview of visio.lign® media



**crea.lign®**  
veneering composite  
A presentation of this  
high-end veneering  
material, its multiple  
applications and unique  
properties.

REF 000577GB



REF 000590GB

**bre.Lux**  
PowerUnit 2  
The LED light-curing  
device is optimised for  
processing crea.lign® - a  
short summary of the  
most relevant informa-  
tion.



REF 000305GB

**neo.lign®**  
Denture teeth  
Design and function of  
the neo.lign® denture  
teeth will fascinate you.  
Discover the variety of  
teeth and select your  
favourite design.

**crea.lign®** - veneering composite


The light-curing composite features impressive natural opalescence and unparalleled light transmission for natural beauty.

**novo.lign®** - veneers

The novo.lign® veneers feature an extended neck and a natural depth thanks to the morphological layer structure. The gold standard in implant prosthetics.

**visio.CAM** - composite blanks

The material benefits from the visio.lign® system are also available for CAD/CAM processing. The blanks from the visio.CAM line are compatible with the entire visio.lign® system.

**neo.lign®** - denture teeth

The neo.lign® denture teeth feature the same shades, designs and materials as the novo.lign® veneers – optimised for fixed/removable restorations.

**bond.lign®** - primer/bonder

The primer/bonder system for a permanent and secure chemical adhesive bond between all system components and common framework materials – including BioHPP® (PEEK).

**visio.paint** - stains

The visio.paint stains enable fast and simple individualisation of veneers.



## visio.lign®

## Info folder

The most relevant materials for the visio.lign® system in a handy info folder.

REF 000274GB



REF 009504GB

## visio.lign®

## Presentation

## of the system

A compact summary of the variety of products of the visio.lign® system.



REF 000234GB

## visio.lign®

## Manual

Step-by-step instructions of the various processing techniques and possibilities of the visio.lign® system.

**crea.lign®**

Veneering composite



create natural beauty  
Red-White Aesthetics

bredent <sup>group</sup>

000651GB-20170124  
Subject to errors and changes

