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(EN) Instructions for use
Tray material UV

Please read these instructions for use carefully before using the product.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Product description

Tray material UV is an acrylate-based light-curing material. Tray material UV is available as bands or as prefabricated plates for the upper jaw.

2. Indication

Tray material UV is suitable for the fabrication of individual trays, functional trays and base plates for setting teeth on.

2.1 Target patient groups

Adults requiring dentures

2.2 Intended users

The product must only be used by dentists, dental technicians and adequately trained dental staff.

3. Contraindication

Tray material UV must not be used intraorally.

If intolerance is suspected, the product must not be used until allergy tests that confirm the absence of an allergy have been completed.

4. Safety recommendations and hazard warnings

Based on the data currently available, the product is not a hazardous substance as defined by EU directives. Nevertheless, the standard precautions for handling chemicals must be taken.

Hazard statements:

EUH204	Contains isocyanates. May produce an allergic reaction.
EUH208	Contains (1-methyl-1,2-ethanediyl)bis[oxy(methyl-2,1-ethanediyl)] diacrylate. May produce an allergic reaction.
EUH210	Safety data sheet available on request.

5. Storage and shelf life information

Storage: Tray material UV must be stored in a tightly closed container at a temperature between 4 °C and 25 °C. **Shelf life:** Under the specified storage conditions, tray material UV can be stored for 5 years.

6. Processing/handling

6.1 Preparation of the object

Insulate the plaster model with Isoplast (REF 54001019) or another plaster-resin insulating agent. Block out undercuts and teeth with Transblock (REF 54001149) or wax.

6.2 Processing of prefabricated plates

Place the plate on the model and shape accordingly.

Please note: do not spread the resin too thinly!

Remove excess resin with a knife or a modelling instrument. Shape the handle from the excess and model it onto the tray plate.

6.3 Processing of band material

Cut off the required amount of tray material UV, place it on the model and shape accordingly. Continue processing in the same way as the prefabricated plates.

Please note: After the required material has been removed from the container, immediately close the container again and protect the remaining material from light.

6.4 Polymerisation

Polymerisation is done occlusally to basally in a suitable light curing unit in the following times:

bre.Lux Power Unit 2	bredent	90 s each
bre.Lux Power Unit	bredent	3 min each
Visio Beta Vario	3M Espe	3 min each
Dentacolor XS	Kulzer	3 min each
Individuo Light Box	Voco	3 min each
Labolight LV-1	Girrbach	3 min each
Spectra SL400	Schütz	3 min each
Spectramat	Ivoclar	3 min each

6.5 Finishing

If additional material has to be applied, this should be done immediately after light curing and prior to finishing since the dispersion layer ensures adequate bonding. Finishing is carried out in the usual way using diamond grinding tools (diacrylic grinding tools) or tungsten carbide burs. If required, the dispersion layer can be removed with alcohol.

7. Technical data

Composition	
Acrylate	~ 24%
Filler	~ 75%
Initiators, stabilisers, pigments	< 1%

Properties based on DIN EN ISO 20795-1 (typical values)	
Flexural strength	50 - 55 MPa
E-modulus	5,000 - 6,500 MPa

8. Other information

The product must only be used by dentists, dental technicians and adequately trained dental staff. For product processing, the instructions for use must be handed to the user as an accompanying document for the insertion or further processing. Original tools and components must be used for processing.

These instructions for use correspond to the current state of the art and our own experience. The product must only be used in the indication described in section 2. The user is responsible for the use of the product. No liability is accepted for incorrect results, as the manufacturer has no influence on the processing of the product. Any claims for damages that are made nonetheless may only be made up to the value of our products.

The user and/or the patient is/are obliged to report any serious incidents occurring in connection with the product to the manufacturer and/or the responsible authority of the Member State in which the user and/or the patient resides.

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