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000197EX -20150313

Stumpflack lichthärtend

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Instructions for use

Die varnish light-curing

Before using the product, please read these instructions for use carefully!

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

The die varnish light-curing is a one-component opaque liquid. This die varnish features high abrasion resistance and therefore withstands high mechanical stress. As a result, the model die is protected against damage when fitting crowns and bridges. The light-curing die varnish is available in blue, yellow, green, red and transparent.

2. Indication

Covering of model dies to protect against damage when fitting crowns and bridges and to create the adhesive gap.

3. Contraindication

Based on the current state of knowledge there are no known contraindications.

4. Hazard and safety information

Signal word:	
Danger	
Hazard statements:	
H225	Highly flammable liquid and vapour.
H316f	Suspected of damaging fertility.
H315	Causes skin irritation.
H319	Causes serious eye irritation.
H317	May cause an allergic skin reaction.
H335	May cause respiratory irritation.
H412	Harmful to aquatic life with long lasting effects.
Precautionary statements:	
Prevention:	
P201	Obtain special instructions before use.
P210	Keep away heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P233	Keep container tightly closed.
P261	Avoid breathing vapours.
P280	Wear protective gloves/eye protection.
Response:	
P303+P361+P353	IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P308+P313	If exposed or concerned: Get medical advice/attention.

Product contains: methyl methacrylate, diphenyl(2,4,6-trimethylbenzoyl)-phosphineoxide, 2-(Dimethylamino)ethylbenzoate

5. Storage and durability

Storage:	
P403+P235	Store in a well-ventilated place. Keep cool.

6. Processing

6.1 Note

To create a cement gap, the die varnish must be polymerized immediately after the application. To harden preparation margins, apply a thin layer of die varnish; remove the excess die varnish with a paper towel before polymerization is started. The application of a thin layer will be sufficient to harden the surface of the die.

The light-curing die varnish allows the individual application of layers and - even on rough undersurfaces - forms a glass-like surface which is resistant to steam-jet blasting and can not be removed from the plaster die.

6.2 Preparation of the object

After the fabrication of the saw-cut model and the plaster dies clean them with compressed air to obtain a dust-free surface. Then apply the desired layer thickness of die varnish light-curing.

Please note: Plaster dies do not need to be pretreated with die hardening and insulating agent.

6.3 Processing

Please note: Bottle must be sealed immediately after use.

Drip die varnish light-curing translucent directly from the bottle onto a mixing block (mixing blocks, REF 330 0114 4) and then apply an evenly thin (approx. 10 µ) layer of die varnish onto the model die (see item 6.5 Polymerization) using a disposable brush, REF 330 0114 2 (brush holder, straight, REF 330 0114 9; brush holder, bent, REF 330 0114 1). If die varnish is applied too richly, it may run off and result in insufficient hardening.

6.4 Layer thickness

The layer thickness can be adjusted by applying the varnish repeatedly and by an intermediate polymerization process of the individual layers. Do not touch the varnish surface between the individual layers. The transparent quality of the die varnish light-curing translucent ensures that markings on the die varnish can be recognized after the application of the varnish. Additionally, the transparent quality allows that uneven application of layers can be recognized easily. Uneven application of layers must be corrected by applying a further layer so that the colour becomes more intense and the layer thickness can be checked. After polymerizing, work can be immediately continued on the varnished die.

6.5 Polymerization

Please note: During light-curing (polymerization) it must be ensured that the varnished surface is as close to the light outlet of the polymerization unit as possible. Dark spots must be avoided. Depending on the respective polymerization unit, the wetted surface must be cured as follows:

6.6 Polymerization times

Dentacolor	approx. 90 - 180 sec
Polylux	approx. 180 - 240 sec
bre.Lux	approx. 180 - 240 sec
Schütz PLC 5000	approx. 180 - 240 sec

7. Symbols

7.1 Symbols according to DIN EN ISO 15223-1

	Reference number
	Attention!
	Manufacturer

7.2 Symbols according to GHS regulation

	GHS02, flame
	GHS 07, Exclamation mark
	GHS08, Health hazard

7.3 Symbols according to NFPA (USA)

	Hazard diamond
--	----------------

8. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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