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Ropak UV

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(EN) Instructions for use

Ropak UV

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

Ropak UV is a light curing resin-colored opaque to cover CoCr objects which are to be veneered with cold-polymerizing, pink resin. The special color of Ropak UV ensures perfect covering of all metal parts.

2. Indication

Opaque for covering CoCr objects (structures) which are to be veneered with pink resin.

3. Contraindication

The product should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.

4. Safety recommendations and hazard warnings

Ropak UV liquid

Signal word:

Danger

Hazard statements:

H225	Highly flammable liquid and vapour.
H361f	Suspected of damaging fertility.
H315	Causes skin irritation.
H319	Causes serious eye irritation.
H317	May cause an allergic skin reaction.
H335	May cause respiratory irritation.
H412	Harmful to aquatic life with long lasting effects.

Precautionary statements:

Prevention:

P201	Obtain special instructions before use.
P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P233	Keep container tightly closed.
P261	Avoid breathing vapours.
P280	Wear protective gloves/eye protection.
P281	Use personal protective equipment as required.

Response:

P303+P361+P353	IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
----------------	--

P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
----------------	--

P308+P313	IF exposed or concerned: Get medical advice/attention.
P362+P364	Take off contaminated clothing and wash before reuse.

Product contains: methyl methacrylate, diphenyl(2,4,6-trimethylbenzoyl)-phosphineoxide, dipentaerythritol pentaacrylate, 2-(dimethylamino)ethyl benzoate

Ropak UV powder

Signal word:

Warning

Hazard statements:

H317	May cause an allergic skin reaction.
H412	Harmful to aquatic organisms, with long lasting effects.

Precautionary statements:

Prevention:

P261	Avoid breathing dust.
P280	Wear protective gloves/eye protection.

Response:

P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
P363	Wash contaminated clothing before reuse.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.

Disposal:

P501	Dispose of contents/container in accordance with local/regional/national/international regulations.
------	---

Product contains: dibenzoyl peroxide

5. Storage and durability information

Storage:

P403+P233	Store container tightly closed in a well-ventilated place.
P403+P235	Store in a well-ventilated place. Keep cool.

Durability: If stored between 4 °C and 25 °C, Ropak UV powder and liquid have an unlimited shelf life.

6. Processing

6.1 Preparation of the object

After casting, finishing and polishing of the chrome cobalt frame, blast the area that is to be covered with Ropak UV with aluminium oxide, grain size 110 µm (units that feature abrasive material circulation are not suitable due to the contamination of the abrasive material since they contaminate the metal surface). After blasting, the metal surface must be cleaned with oil-free compressed air or a steam-jet unit.

6.2 Mixing

Mix Ropak UV powder and liquid (e.g. on a mixing block) to obtain a liquid consistency.

Please note: Do not mix Ropak UV so that the consistency is too thick since this will result in incomplete curing.

6.3 Applying

Apply Ropak UV in two thin, covering layers using a disposable brush and carry out intermediate polymerization. Then polymerize in the light curing unit as described under item 6.4.

6.4 Technical data

Denta-Color XS	90 sec.
Schütz PLC 5000	5 min.
bredent PolyLux	5 min.
bredent bre.Lux PowerUnit/ bre.Lux PowerUnit 2	180 sec. intermediate polymerization, 360 sec. end polymerization

Please note: The object should always be placed as close as possible to the light source. The luminous intensity, and thus the effect decreases considerably if the distance to the light source increases, so that Insufficient polymerization may result. Shade spots must be avoided. Ropak UV hardens between 470 to 520 nm.

7. Sources of errors

7.1 Formation of „clouds“ and discoloration in the resin

Error:	Ropak UV has not been cured long enough in the light-curing unit.
	The layer of Ropak UV that has been applied is too thick.
	The material has been cured in an unsuitable light-curing unit.

7.2 Despite sufficiently long curing times Ropak UV still exhibits greasiness on the surface.

Error: The layer of Ropak UV is too thick.

Solution: Remove Ropak UV by sandblasting the frame (object); mix it to obtain a thinner consistency and apply anew.

8. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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