# Brückenteilungsgeschiebe OC

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## Instructions for use

# OC Bridge Sectioning Attachment

## 1. Description of the product and range of applications

Cylindrical bridge sectioning attachment - screw-retained through the occlusal surface.

## 2. Processing

### 2.1 Caution

**Please note:** At least 4.3 mm must be available vertically for the bridge sectioning attachment, including screw and locking ring. The attachment and screw cannot be shortened more than this. If less space is available, we recommend using the „Custom Bridge Sectioning Attachment“ (REF 430 0735 0) and the „Tool set“ for custom threaded retainers (REF 330 0060 0). This enables the screw to be set at the required angle to the attachment stud.

### 2.2 Preparation

Wax up the crown for the bridge sectioning attachment using standard methods. Use the paralleling mandrel (REF 360 0115 7) to position the bridge sectioning attachment on the crown, according to the angle of insertion of the second bridge section. The underside of the bridge sectioning attachment can be adapted to the papillae.

**Caution:** The plastic section can only be reduced to the threaded sleeve, i.e. 3.3 mm.

### 2.3 Prior to investing

Screw in the retention rod to hold the threaded sleeve in the correct position in the casting ring. Before screwing the rod into the sleeve, coat it with a thin layer of colloidal graphite (Order No.: 540 0070 6). This prevents the rod burning into the investment and enables it to be removed easily after casting. When screwing in the retention rod, ensure that it is tightened to the stop.

### 2.4 Investing and casting

The cast-on threaded sleeve of the bridge sectioning attachment consists of a high gold content (Au, Pd, Pt) alloy (melting range: 1320°C - 1460°C). Therefore, the liquidus of the alloy being cast on must not exceed 1320°C. Invest and cast using standard procedures and according to the investment material, alloy and casting machine manufacturers' instructions.

### 2.5 Devesting

After casting, devest and clean as usual. Do not remove the retention rod until the castings have been sandblasted.

**Caution:** The retention rod can only be used once as it will not function properly for a second restoration.

### 2.6 Trimming and milling

Trim the crown and attachment as usual. The section connecting the attachment to the coping is rounded and has a diameter of 1.0 mm. Therefore, the attachment can be trimmed with an ISO size 010 cylindrical cutter if necessary. The restoration should be polished according to the alloy used.

### 2.7. Fabricating the outer section

Screw the locking ring into place with the titanium screw. Grind in both components according to the bite. The ridge around the locking ring marks the level to which the locking ring and titanium screw can be reduced. Neither should be reduced further as, otherwise, the hexagonal socket in the screw will no longer be sufficient to tighten the screw with the screwdriver. Once the locking ring has been reduced, the hexagonal socket must remain 0.6 mm deep. Then coat the screw, locking ring and bridge sectioning attachment with fully combustible resin (e.g. Pi-Ku-Plast, REF 540 0017-) and wax up the outer section. Before continuing waxing up the bridge, raise the outer section and reposition it.

### 2.8 Investing and casting

The cast-on locking ring of the bridge sectioning attachment consists of a high gold-content alloy (melting range: 1320°C - 1460°C). Therefore, the liquidus temperature of the alloy being cast against them not exceed 1320°C. Invest and cast using standard procedures and according to the investment material, alloy and casting machine manufacturers' instructions.

## 3. Safety hints

### 3.1 Protective clothing

Always wear protective clothing when working.

### 3.2 Working with Pi-Ku-Plast

When processing Pi-Ku-Plast, the procedure described in the instructions for use must be observed.

### 4. Please note

The information contained in these instructions for use is always updated according to the latest knowledge and experience. Therefore we recommend to read the instructions for use again before using a new package. These non-binding user recommendations are based on our own experience. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since bredent does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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