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Plaster-resin insulating agent Isoplast ip
- Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.
- 1. Description of the product**

Isoplast ip is an alginate-based insulating agent. It ensures dependable separation of plaster against resin. Isoplast ip is available in a 750 ml bottle.

**2. Indication**

Isoplast serves to separate plaster models against hot- and cold-curing resins as well as silicones. In particular the indication range includes the fabrication of full dentures in the casting, condensing and injection technique and partial dentures.

**3. Contraindication**

Based on the current state of knowledge there are no known contraindications.

**4. Safety recommendations and hazard warnings**

**Signal word:**

Attention

**Hazard statements:**

H317 May cause an allergic skin reaction.

**Precautionary statements:**

**Prevention:**

P261 Avoid breathing/vapours.  
P280 Wear protective gloves/eye protection.

**Response:**

P302+P352 IF ON SKIN: Wash with plenty of water.  
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.  
P363 Wash contaminated clothing before reuse.

**Product contains:** mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2-Hisothiazol-3-one [EC no. 220-239-6] (3:1)

**5. Storage and durability**

**Disposal:**

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

**6. Processing**

For processing, Isoplast ip is filled into a small bottle that can be closed and then it is spread out of the bottle onto the plaster models to be separated using a clean brush. Thus contaminating of Isoplast ip in the storage bottle is avoided. In order to obtain a homogeneous insulating layer, any wax residues must be removed from the model and the model must be soaked. Apply Isoplast ip onto the warm or lukewarm plaster using a soft brush. Puddles and uneven layer thicknesses can be easily recognized due to the red color. Excess insulating agent is absorbed with a clean brush. If a single layer of the insulating agent is applied, the layer thickness is approx. 200 µm. Isoplast ip has dried after approx. 6 minutes and seals the surface of the model. Then the plaster model can be processed further.

**7. Additional information**

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products. Enamel that has already been etched must not be contaminated. If it is contaminated, etching, rinsing with water and drying must be repeated as described above. The etching gel must not come in contact with skin, mucous membranes or eyes. In case of accidental contact, rinse with copious amounts of water.

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bredent GmbH & Co. KG