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Exakto-Form Modellkunststoff

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(EN) Instructions for use

Exakto-Form model resin

Before using the product, please read these instructions for use carefully!

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

Exakto-Form model resin is a two-component material based on polyurethane (PU). The extremely high fluidity of Exakto-Form model resin allows pouring of impressions with high precision and without the formation of bubbles. True-to-the-original and accurate reproduction of surfaces and perfect dimensional stability as well as superior precision of fit without expansion are distinctive features of this model resin.

2. Indication

Exakto-Form model resin is suitable for the fabrication of sawcut models and master models and for individual die segments.

3. Contraindication

Exakto-Form model resin must not be used for hydrocolloid impressions. The product should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.

4. Safety recommendations and hazard warning

Exakto-Form insulating liquid

Signal word:

Danger

Hazard statements:

H226	Flammable liquid and vapour.
H304	May be fatal if swallowed and enters airways.
H319	Causes serious eye irritation.

Precautionary statements:

Prevention:

P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P233	Keep container tightly closed.
P280	Wear protective gloves/eye protection.

Response:

P301+P310	IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician.
P303+P361+P353	IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
----------------	--

P331 DO NOT induce vomiting.

Product contains: Alkanes, C9-12 iso

Exakto-Form model resin Component A

Signal word:

Warning

Hazard statements:

H319	Causes serious eye irritation.
H411	Toxic to aquatic life with long lasting effects.

Precautionary statements:

Prevention:

P273	Avoid release to the environment.
P280	Wear protective gloves/protective clothing/eye protection/face protection.

Response:

P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
	If eye irritation persists: Get medical advice/attention.

P337+P313	
Disposal:	
P501	Dispose of contents/container in accordance with local/regional/national/international regulations.

Exakto-Form model resin Component B

Signal word:

Danger

Hazard statements:

H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335	May cause respiratory irritation.
H351	Suspected of causing cancer.
H373	May cause damage to organs through prolonged or repeated exposure.
H411	Toxic to aqastic life with long lasting effects.

Precautionary statements:

Prevention:

P260	Do not breathe dust/fume/gas/mist/vapours/spray.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P284	[In case of inadequate ventilation] wear respiratory protection.

Response:

P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P304+P341	IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P342+P311	If experiencing respiratory symptoms: Call a POISON CENTER/doctor.

Disposal:

P501	Dispose of contents/container in accordance with local/regional/national/international regulations.
	Contains isocyanates. May produce an allergic reaction.

Product contains: diphenylmethanediisocyanate,isomeres and homologues, Bis(isopropyl)naphthalin

5. Storage and durability

5.1 Storage

P403+P235 Store in a well-ventilated place. Keep cool.

5.2 Durability

After receiving Exakto-Form model resin, use up the material quickly. In case of prolonged storage, the components may thicken which, however, does not cause a loss in quality. In such cases, stir the components thoroughly and completely.

6. Disposal

Disposal:

P501	Dispose of contents/container in accordance with local/regional/national/international regulations.
------	---

If it is not processed (original condition), Exakto-Form model resin must not reach the sewage systems/bodies of water and the soil. Hardened Exakto-Form model resin can be disposed of together with the domestic waste.

7. Processing

Please note: Gloves must be worn when processing Exakto-Form model resin.

7.1 Mixing the components

Each component must be stirred separately with a stirring stick (REF 39000310) prior to mixing to obtain a homogeneous mixture. Mix the sediment completely. Add component B (resin) into component A; empty the tin completely. Mix the components thoroughly for approx. 30 seconds until a uniform color is obtained. Pour the mixed components swiftly into the impression. The processing time span is approx. 2 to 3 minutes. Exakto-Form model resin can be removed 30 minutes after mixing. Final hardness is achieved after 90 minutes.

To exclude premature reaction of component B after opening (upon the penetration of air moisture), processing should be carried out in a swift manner.

7.2 Fabrication of an arch

Exakto-Form model resin is suitable for silicone and polyether impressions. Clean the impression and dry it thoroughly. Seal the dorsal open area with a wax sheet to avoid any leakage of Exakto-Form model resin when pouring the impression. Spray a thin coat of Exakto-From insulating liquid onto polyether impressions to avoid adhesion. Mix Exakto-Form model resin as described under 7.1 and pour it slowly (fine stream) into the impression. Wait at least 30 minutes to remove the arch and continue processing in the usual way after 90 minutes.

Exakto-Form model resin is processed like plaster using tungsten carbide burs, a pin drill unit and diamond separating discs. Use Exakto-Form insulating liquid to separate the arch and the model base.

7.3 Fabrication of a master model

Exakto-Form model resin is suitable for silicone and polyether impressions. Clean the impression and dry it thoroughly. Apply functional margin wax (REF 43001500) around the vestibular area of the impression to obtain a border and functional margin. Spray a thin coat of Exakto-Form insulating liquid onto the base former to avoid adhesion. Pour the mixed components (see 7.1) thinly into the impression and the base former. The impression should only be placed on the base former when the Exakto-Form resin has reached a paste-like consistency (risk of leaking!).

7.4 Fabrication of a die segment

Exakto-Form model resin is suitable for silicone and polyether impressions. Clean the impression and dry it thoroughly. Apply wax to the border of the corresponding segment. Spray a thin coat of Exakto-Form insulating liquid onto polyether impressions to avoid adhesion.

The components are mixed separately and then each component is filled separately into the dosing/measuring syringes (REF 39000360). Add equal portions of the desired quantity into the mixing cup (REF 39000300) and mix with the stirring stick. Then fill the material into the segment in the impression. The syringes must be closed airtight and the content must be processed within 2 days.

8. Technical data

Mixing ratio	1:1
Processing and hardening temperature	18 to 23 °C
Processing time span	2 to 3 minutes
Time until removal of the impression	30 minutes
Time to achieve final hardness	90 minutes
Density	1.6 kg/dm³
Shore hardness D (DIN 53505)	80
Flexural strength (DIN 53452)	50 N/mm²
e-modulus (DIN 53452)	3900 N/mm²
Linear shrinkage	0.17 %

9. Additional information

The product may only be used by dentists, dental technicians and adequately trained dental staff. Only original tools and components may be used for processing. These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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