Master-Pin-System

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Instructions for use

Master-Pin-System

Before using the product, please read these instructions for carfellyu use.

1. Indication range

The Master-Pin system, Master-Pin (brass) and Master-Pin sleeve (plastic), is used for exact repositioning of sawed individual dies on the model base of the sawcut model.

2. Processing

2.1 Preparation of the arch

The arch into which the holes for the pins are to be drilled must be ground flat. A plaster trimmer with diamond trimming wheel is perfectly suitable for this purpose.

The trimmed base of the arch can be optimized with wet sanding paper (e.g. grain size 600). Let the arch dry well to avoid smearing of the plaster during milling and drilling.

Prior to drilling and placing the Master-Pins in the arch, grind the the lingual resp. palatal inner side of the arch slightly conical (6°) towards the base. This will simplify subsequent removal of the arch from the model base and also serves to firmly position the individual dies on the model base and protect them against torsion.

2.2 Drills

Depending on the pin drilling unit the diamater of the drill hole may vary due to imprecisions of the collect chuck. Therefore three different pin drills are available for the Master-Pin system.

2.2.1 Master-Pin Diatit tungsten carbide step drill standard/green ø 1.5/2 mm

The Master-Pin assortment includes the Master-Pin Diatit tungsten carbide step drill standard/green ø 1.5/2. In most cases, the Master-Pin Diatit tungsten carbide step drill ø 1.5/2 will drill precise drill holes to receive the Master-Pins.

2.2.2 Master-Pin Diatit tungsten carbide drill special/yellow ø 1.5/2 mm

If it is difficult to insert the pins into the arch due to the concentricity tolerance of the pin drilling unit, the Master-Pin Diatit tungsten carbide step drill special/yellow ø 1.5/2 which features a larger diameter (+ 0.01 mm) can be used. This drill will prepare a larger drill hole so that the Master-Pins can be inserted and glued in more easily.

2.2.3 Master-Pin Diatit tungsten carbide step drill special/red ø 1.5/2 mm

If the drill holes in the arch are too large due to the concentricity tolerance of the pin drilling unit with the Master-Pin Diatit tungsten carbide step drill standard/green ø 1.5/2 – the Master-Pin does not have a stable, fixed position in the drill hole – the Master-Pin Diatit tungsten carbide step drill special/red ø 1.5/2 which features a smaller diameter (- 0.01 mm) can be used. This drill will prepare a smaller drill hole.

2.3 Drilling the pin holes

Precondition for this work is the correct adjustment of the Master-Pin Diatit tungsten carbide step drill standard / green ø 1.5/2 mm (Order No. 360 0119 2) in the pin drill unit.

Recommendation: Prior to drilling, place dry arches shortly on a wet cloth with the basal side facing downward. Then blow them dry and drill the drill holes (correct position must be ensured) immediately using the pin drilling! unit. Breaking out of plaster around the drill holes will be avoided – precise drill holes are obtained.

The boundary line for drilling at the Master-Pin must be flush with the surface of the arch to receive the exact drilling depth. We recommend to prepare two drill holes for each die:

1. Drill hole = central fissure; 2. Drill hole = at a distance of 2 to 3 mm toward the lingual resp. palatal area.

2.4 Glueing in the Master-Pin

Add a drop of cyanoacrylate adhesive (super glue) onto a suitable object (e.g. wax plate). Only the surface of the face of the glueing shaft of the Master-Pin is dipped in the cyanoacrylate adhesive. The Master-Pin is inserted up to the stop into the drill hole at a right angle towards the base of the arch.

2.5 Separating

After glueing in the Master-Pins, apply Master-Sep (Order No. 520 0029 0) plaster/plaster separating liquid onto the die using a brush.

2.6 Sliding on the Master-Pin sleeves

Slide the Master-Pin sleeves onto the Master-Pins so that they are seated flush with the base of the arch. Should two master pins be positioned very closely to each other, the two Master-Pin sleeves are slid on with the flattened sides facing each other.

2.7 Preparing the base

Prepare several arches using liquid plaster and e.g. the Master-Split system (Order No. 360 0124 M) without the formation of bubbles.

Since the Master-Pin sleeves are 0.5 mm higher than the Master-Pins, a uniform, defined height is achieved. All Master-Pin sleeves can be seen at the model base.

2.8 Trimming

Finally, trim the sides of the model to ensure that the arch can be easily removed from the model base.

3. Safety recommendations and hazard warnings

When working with the pin drilling unit, the safety recommendations and hazard warnings of the respective manufacturer must be observed; the same applies to handling of the cyanoacrylate adhesive (super glue).

4. Storage and durability

The Master-Pin system must be stored at room temperature under dry conditions.

5. Other information

The information contained in these instructions for use is always updated according to the latest knowledge and experience. Therefore we recommend to read the instructions for use again before using a new package. These instructions for use refer to Master-Pin system, Order No. 360 0122 6, 360 0122 5 and 360 0119 2.

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