Gebruiksaanweisung

Kugelfixator

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Instructions for use

Stud fixator

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

The stud fixator consists of a ceramic stud and a cavity-filling silicone that acts like a buffer. When fabricating new telescopic restorations or when repairing them, the stud fixator can be used to increase or restore the friction.

2. Indication

The stud fixator is a temporary tool to restore lost friction of telescopic and conical crowns and for new telescopic and conical crowns. A sufficient coping thickness (0.4 mm) of the primary construction is essential for processing.

3. Contraindication

Based on the current state of knowledge there are no known contraindications.

4. Safety recommendations and hazard warnings

Please refer to the safety data sheet for safety recommendations concerning the product. The lot numbers of all components in use must be documented to ensure traceability and to handle complaints.

bredent

GmbH & Co.KG · Weissenhorner Str. 2 · 89250 Senden · Germany
 Tel. (+49) 0 73 09 / 8 72-22 · Fax (+49) 0 73 09 / 8 72-24
 www.bredent.com · e-mail info@bredent.com

5. Storage and durability

Storage No special storage conditions are required for the stud fixator.

Durability

There are no durability limits for the stud fixator.

6. Processing

6.1 Processing in the laboratory

A removable die, preferably made of resin, and a working model are fabricated in the laboratory. Then a matrix of the resin saddle is produced to ensure the accurate position of the resin saddle when repositioning later on. The resin saddle is removed. A coping produced from thermoforming foil is placed on the resin die; it serves as drilling template. The precise position for the integration of the stud fixator is determined at the secondary crown and a hole for the stud fixator is drilled into the secondary crown (Ø 2.4 mm) using a tungsten carbide bur. The restoration is placed on the model again and the accurate position of the recess in the resin die becomes evident by slight manual turning of the Rapidly 1.4, (REF H001 NH 14), through the drill hole in the secondary element. This recess must have a max. depth of 0.3 mm in the resin die. The coping produced from the thermoforming foil is placed on the die and the drilling position is marked. A hole is drilled at the marked position that serves as drilling template for the dentist later on. The stud fixator is fixed in the secondary crown and attached to the CoCr structure using DTK adhesive, (REF 540 0010 6). The resin saddle is fastened again with the help of the matrix and polished.

6.2 Processing in the mouth

The dentist places the coping produced from the thermoforming foil onto the corresponding primary element and drills a groove with a maximum depth of 0.3 mm using the Rapidly bur 1.4, (REF H001 NH 14). The denture can be integrated.

7. Technical data

Dimensions: 2.4 x 3.4 mm
 Material: titanium, silicone, ceramic

8. Symbols

LOT

Lot number

9. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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