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compoForm UV

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Instructions for use

compoForm UV

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

compoForm UV is a light-curing dental composite (for modeling) based on methacrylates. The product can be directly applied from the syringe and hardened using standard manual polymerization lamps (UV and LED). compoForm UV can be combined with all modeling waxes and burns out without leaving any residue.

2. Indication

- Modeling post restorations (extraoral)
- Modeling casting elements in the milling and attachment technique
- Extension of CoCr restorations
- Interlocking (connecting) bridges for soldering
- Blocking out undercuts when preparing models

3. Contraindication

compoForm UV must not be used in the oral cavity.

4. Safety recommendations and hazard warnings

Signal word:
not applicable

Hazard statements:

H412 Harmful to aquatic life with long lasting effects.

Precautionary statements:

Prevention:

P273 Avoid release to the environment

Disposal:

P501 Dispose of contents/container in accordance with local / regional / national / international regulations.

5. Storage and durability

compoForm UV must be stored under dry conditions. Exposure to direct sunlight and heat must be avoided. Please observe the expiration date printed on the label.

6. Properties

compoForm UV is directly applied from the syringe. The material features a readily flowing and yet stable consistency. Thanks to direct hardening during the application, larger models can also be prepared easily. A perfect casting result is ensured thanks to the fact that compoForm UV burns out without leaving any residue.

7. Processing

After the cap has been removed, one of the enclosed cannulas is screwed on to ensure highly precise dos-ing when the material is applied. The cannulas are intended for single use only. compoForm UV can now be applied in thin layers. To avoid the formation of stress in the model, the layers that are applied should not be thicker than 2 mm and hardened with the polymerization lamp immediately during the application. To ensure perfect bonding between the layers, the inhibition layer that is formed during the individual polymerization steps must not be removed. If compoForm UV is to be applied subsequently, the model must be free from dust and grease. compoForm UV can be finished with standard instruments/tools.

8. Preheating of the furnace

Recommended preheating for the use of compoForm UV:

8° C/min.	→	270° C
Holding time	=	30 minutes
8° C/Min.	→	580° C
Holding time	=	30 minutes
8° C/Min.	→	Final temperature
Holding time	=	60 minutes

9. Symbols

9.1 Symbols according to DIN EN ISO 15223

	Reference number
	Lot number
	Can be used until YYYY-MM
	Observe accompanying documents
	Protect against sunlight
	Temperature limits; upper limit: 25° C, lower limit: 4° C
	Manufacturer

9.2 Symbols according to NFPA (USA)

	Hazard diamond
--	----------------

10. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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