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bre.Lux Power Unit, bre.Lux LED N, bre.Lux Flexholder, bre.Lux Updown

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Instructions for use

bre.Lux Power Unit, bre.Lux LED N, bre.Lux Flexholder, bre.Lux UpDown

Before using the product, please read these instructions for use carefully!

1. Application and use

bre.Lux Power Unit is a light-curing (polymerization) unit for the dental laboratory or the dentist practice and is used for curing light- and / or dual-curing pastes and liquids. Light emitting diodes (LEDs) with different wavelengths are used for curing and hence the bre.Lux Power Unit covers the wavelength range of 370 to 500 nm and is therefore suitable for all standard dental materials. The 21 power LEDs of the unit have a service life of 20,000 hours without any significant change to the performance. Moreover bre.Lux Power Unit is suitable to accommodate a hand lamp and to supply power to the hand lamp which allows the use of both hands if a flexible hose is attached.

2. Product types

2.1 bre.Lux Power Unit

The scope of delivery include the bre.Lux Power unit polymerization unit, a bre.Lux LED N hand lamp and the respective coiled cord, a flexible hose, a power cord and a bre.Lux UpDown object table. Please make sure that all components are contained in the package.

Attach the flexible hose (1) to the rear (on the top left side next to the main switch (2)). Use the coiled cord (3) to connect the hand lamp to the unit (rear of the unit – top right (4)). Once the flexible hose has been aligned, the hand lamp can be fixed in the receiving ring.


The unit must have a minimum distance of at least 10 cm to the rear wall to ensure adequate circulation of air. The bre. Lux UpDown object table has 2 different foot lengths to be able to reduce/increase the distance of the objects to be polymerized to the LEDs mounted in the top section of the chamber.

Unit with flexible hose

Weight: approx. 3200 g

Dimensions: 300 x 225 x 105 mm

Internal dimension: Polymerization chamber 180 x 110 x 75 mm

2.2 Description of keys and functions of keys

- (1) Display
- (2) 20 seconds with 50% power
- (3) 40 seconds with 50% power
- (4) Permanent reduction to 70% power
- (5) 90 seconds – full power
- (6) 120 seconds – full power
- (7) 180 seconds – full power
- (8) 360 seconds – full power
- (9) Continuous operation at full power

2.3 Keys of hand lamp

- (1) Start / Stop + continuous function
- (2) 15 seconds
- (3) 30 seconds

3. Putting into operation

bre.Lux Power Unit

- ➔ Well ventilated space, distance of at least 15 cm to the rear panel.
 - ➔ First connect the power cord with the unit and then to AC power source.
 - ➔ Use the main switch (rear of the housing) to switch the unit on.
 - ➔ The display shows „bre.Lux Power unit“.
 - ➔ The fan will start when the unit is put into operation regardless of any key that has been previously selected. No air turbulence will be formed inside the chamber!
 - ➔ All default settings are shown on the display. If the start phase is selected the sequence to full power is displayed.
 - ➔ Any key can be pressed to abort all programs.
 - ➔ When the lid is opened, the program will be interrupted and will start at the point of time when it was interrupted
 - ➔ Corrections can be easily made by pressing a different key and confirming it.
- The desktop polymerization unit features a combination of programmed time intervals and continuous operation as well as different power levels which can be combined with one another.

Time interval	Direct selection of the time interval	20 Sec. prepolymerizing with 50 % power	40 Sec. prepolymerizing with 50 % power	Time interval or permanent reduction of power to 70 %
90 Sec.	90 Sec. 100% (5)	20 Sec. 50 %, (2)+(5) 70 Sec. 100 %	40 Sec. 50 % (3)+(5) 50 Sec. 100 %	90 Sec. 70 % (4)+(5)
120 Sec.	120 Sec. 100% (6)	20 Sec. 50%, (2)+(6) 100 Sec. 100 %	40 Sec. 50 % (3)+(6) 80 Sec. 100 %	120 Sec. 70 % (4)+(6)
180 Sec.	180 Sec. 70% (7)	20 Sec. 50 %, (2)+(7)	40 Sec. 50 % (3)+(7)	180 Sec. 70 % (4)+(7)
180 Sec.	180 Sec. 100% (7)	160 Sec. 100 %	140 Sec. 100 %	
360 Sec.	360 Sec. 100% (8)	20 Sec. 50%, (2)+(8) 340 Sec. 100 %	40 Sec. 50 % (3)+(8) 320 Sec. 100 %	360 Sec. 70 % (4)+(8)
Continuous operation	Permanently 100% (9)			Permanently 70% (4)+(9)

If the operating time is directly selected by pressing the respective key (with the desired time interval), polymerization with 100% power is carried out. If prepolymerizing for 20 or 40 sec or continuous operation with reduced power is desired, the keys on the left side can be pressed. Then the keys on the right side are used to select the polymerization time and the unit will start the polymerization process. The countdown is shown on the display.

Generally, LED's will not be switched off / deactivated; all LEDs are continuously activated! It is particularly useful to select the prepolymerization option to cure large material quantities. Moreover this will avoid / reduce contraction towards one direction (e.g. toward the light). Reducing the power to 70 % is recommended for example when using light-curing splint materials or if the material is to be cured from one side first and then from the other in a second step. The prepolymerizing function is also useful for primer and bonder combinations if the reaction time of the primer is to be extended.

4. bre.Lux LED N with coiled cord (included in the scope of delivery)

- ➔ As require, select 15 s or 30 s and press Start key.
- ➔ Continuous operation without selecting a default time, press Start key.
- ➔ Abort any function with the Start / Stop key.
- ➔ The different time functions – 15s, 30 s and continuous operation – blue can be selected by pressing the keys 15 s = green LED and 30 s = red LED twice. The setting selected last remains activated until the next change is performed.

Prepolymerization hand lamp with carefully chosen wavelength range of 370 to 500 nm; power supply using the enclosed coiled cord to be attached to the rear of the housing of the bre.Lux Power Unit.

For curing crea.lign and light-curing die varnishes or for the fixation of veneers using combo.lign.

5. bre.Lux LED N (with plug-in power unit) / REF 140 0097 1

Prepolymerization hand lamp with carefully chosen wavelength range of 370 to nm: for single use, use enclosed plug-in power unit for power supply; other features as described above.

Weight: approx. 225 g

Dimensions: 170 x 25 mm

Note:

To be able to operate the hand lamp included in the scope of delivery of bre.Lux Power Unit (REF 140 0097 0) independent of the basic unit, a separate plug-in power unit (REF 140 0097 9) is required. When using older hand lamps which do not have an integrated coiled in the handle of the hand lamp, a special coiled cord (REF 993 0001 0) must be ordered in addition to the separate plug-in power unit. To enable working with both hands in various positions under the hand lamp and to ensure safe storage of the device when it is not used, we recommend the high-quality Flexholder which ensures safe hold and use thanks to the gooseneck design and facilitates the work process significantly.

6. bre-Lux Flexholder / 140 0097 3

Stainless steel plate with flexible hose (length: 30 cm) and receiving ring to allow working with both hands.


Weight: approx. 560 g

Dimension of round plate: 120 x 6 mm

Length of flexible hose: approx. 230 mm

7. bre.Lux UpDown / REF 140 0097 4

Mirror-polished stainless steel table with leg heights of 20 and 35 mm to increase/reduce the distance of the objects to be polymerized to the LEDs mounted in the upper section of the chamber. The mirror finish ensures excellent reflection of incident light rays. The large number of holes and oval cutouts allow placement of individual dies and bridge structures.

Weight: approx. 48 g

Dimensions: 89 x 89 x 55 x 1 mm

Material: stainless steel

8. Hazard warnings and safety information when operating the LED polymerization units

Do not use a damaged unit if the glass of the LED is missing or defective. The light beam must never be directed towards exposed soft tissue to avoid any injuries or irritations. The light beam must never be directed towards the eyes. Light that is reflected may also cause injuries to the eyes. Use the protective eye shield or safety goggles with a suitable light filter. If the ambient temperature and other sources of heat cause overheating, the unit aborts the program and the program is continued automatically once the unit has cooled down.

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GmbH & Co.KG · Weissenhorner Str. 2 · 89250 Senden · Germany