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Friction Splint FS1

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(EN) Instructions for use

Friction Splint FS1

Before using the product, please read these instructions for use carefully!

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

Friction Splint FS1 is used for partially removable dentures and serves to restore the masticatore function as a connection element for primary and secondary structures. If required, the screw can be exchanged and the partially removable part of the denture can be cleaned.

Friction Splint FS1 consists of a special plastic sleeve and a splint with a thread to ensure safe retention in the patient's mouth. The short screwdriver is used for easy removal of the splint for control and cleaning purposes.

2. Indication

Horizontal screw retention of bridges on residual teeth or implant-supported structures made of the following materials.

- Ceramics
- Zirconium
- High-performance polymers
- Precious metal alloys
- NPM alloys

2.1 Target patient group

The patient group consists of adults with residual teeth in the upper/lower jaw, residual teeth in combination with implants or patients with implant-supported restorations in the upper/lower jaw.

2.2 Intended users

The product may only be used by dentists, dental technicians and adequately trained dental staff.

3. Contraindication

Friction Splint FS1 is contraindicated in the following cases:

- the product should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.

4. Safety and hazard information

The lot numbers of all products used should be recorded to ensure traceability and to handle complaints.

The conditions for proper storage (see item 5) and handling of the products must be observed. Improper storage can influence the product properties and result in failure of the restoration.

When using the product, all parts must be protected against aspiration and accidental swallowing. Aspiration and swallowing of products may cause infections and physical injury.

All non-sterile parts must be disinfected prior to integration into the patient's mouth. All assembled parts need to be disassembled into their individual components.

To avoid cross contamination, the manufacturer's instructions on reprocessing (disinfecting and sterilizing) must be observed.

If abutments are removed for cleaning purposes, new retention screws must be used. Reuse of the retention screws is contraindicated.

The retention screws used in the laboratoy are not suitable for clinical use since they are subject to stress in the laboratory; as a result long-term and reliable fit can no longer be guaranteed.

Friction Splint FS1 screw is intended for single use. However, if this product is used several times, the following risks may arise: cross contamination, malfunction, inaccurate fit, etc.

The specified torque of the abutment systems (item 7.1) must not be exceeded, as overtightening the abutment can lead to screw fracture.

Prevention:

The safety information must be read and understood prior to use. Use/wear required personal protective equipment.

Product contains:

Auxiliary instruments may contain nickel. A list of the materials used can be found in the table under item 7 Technical data.

5. Storage and durability

The product is supplied in the non-sterile condition. Dry and dust-free storage in the original packaging at room temperature is required.

Storage:

Please observe storage information on the label.

Disposal:

Products, instruments and accessories, which do not have any special requirements from the aspect of infection protection as far as collection and disposal are concerned, can be disposed of together with household waste if necessary. Regional, national and international regulations must be observed.

CAUTION

The maximum wearing period of the Friction Splint FS1 sleeve is 12 months. The dentist must inform his patient that he must visit the practice every 12 months to have the splinted structures checked.

6. Processing/ handling

CAUTION

Generally, we recommend thorough surgical and prosthetic planning and selection of the prosthetic components depending on the situation and type of restoration.

All non-sterile parts must be disinfected prior to integration into the patient's mouth. All assembled parts need to be disassembled into their individual components.

The manufacturer's instructions on the use of the disinfectant must be observed.

6.1 Use

6.2 Implant design

A wax-up is prepared to ensure accurate placement of the splinting. The wax-up is fixed with a silicone key to achieve control and orientation for the subsequent working steps. The drill hole should be positioned in a way to ensure that the splint can be easily exchanged in the mouth. If the exact position has been determined, the Diatit-Multidrill Ø 2.0 mm (REF 33000720) is used to drill the hole.

CAUTION

Make sure to use sharp drilling tools, drilling oil and a max. speed of 5000 rpm to prepare a splint hole with an exact diameter of 2.0 mm. Inaccurate splint holes may affect the proper function of the Friction Splint FS1.

When using hard alloys, previous center punching is required; select a drill with a small diameter to start drilling and gradually expand the drill hole to 2.0 mm. After drilling, the auxiliary modelling element (REF 45000083) is placed into the drill hole and the secondary structure is modelled.

After modelling, the auxiliary modelling element is removed and investing is carried out.

Prior to investing, the spacer (REF 45000082) can be placed into the drill hole to facilitate investing and ensure the quality of the splint hole. After casting, devest and sandblast the secondary structure; the spacer is also sandblasted. The secondary structure is fitted and completed.

6.3 Crown design

Splinting of crowns is done as described under § 6.1.

6.4 Fitting in the sleeve and the splint

Sleeve and splint should always have the same length (in accordance with the splint hole).

CAUTION

When shortening the splint and the sleeve, always start from the thread section and never from the head.

The sleeve is inserted into the splint hole together with the splint screw (screwed in up to half of its length) and slight friction is adjusted. The splint screw is screwed in congruently to the sleeve and the splint hole.

6.5 Cleaning and disinfecting

Before inserting the Friction Splint FS1, it is recommended to place the sleeve and the splint into a standard disinfectant. The immersion times given by the manufacturers must be adhered to. Sleeve and splint may only be removed and disinfected by the dentist/surgeon. It is recommended to integrate the patient into a recall system to ensure the maximum wearing period of the sleeve of one year.

6.6 Zirconium

For all-ceramic objects it must be ensured that the splint hole has a diameter of 2.0 mm after the sintering processing and after the ceramic firing processes. The shrinkage factor of the ceramic needs to be taken into account.

7. Technical data

REF:	Product:	Quantity:	Material:	Ø:	Length:
45000080	Friction Splint FS1 sleeve	1 piece	POM	2,0 mm	8 mm
45000081	Friction Splint FS1 pin	1 piece	Titanium grade 5 3.7165	M1.4	8 mm
45000082	Friction Splint FS1 Spacer Ø 2.0 mm	2 pieces	Ceramic	2.0 mm	15 mm
45000083	Friction Splint FS1 Modeling aid Ø 2.0 mm	2 pieces	Stainless steel 1.4305	2.0 mm	15 mm
45000084	Friction Splint FS1 sleeve	10 pieces	POM	M1.4	8 mm
45000085	Friction Splint FS1 pin	10 pieces	Titanium grade 5 3.7165	M1.4	8 mm
45000086	Friction Splint FS1 Spacer Ø 2.0 mm	10 pieces	Ceramic	2.0 mm	15 mm
45000087	Friction Splint FS1 Modeling aid Ø 2.0 mm	10 pieces	Stainless steel 1.4305	2.0 mm	15 mm
45000088	Friction Splint FS1 set 13 Parts	1 set			

8. Products for joint use

CAUTION

Prior to the use of the products for joint use, observe the corresponding instructions for use.

Abutments

8.1 Restrictions with regard to combination

Only original tools and components may be used for processing.

9. Additional information

These instructions for use need to be handed over to the user as an accompanying document for insertion or further processing. Only original tools and components may be used for processing.

CAUTION

When inserting the restoration and during the patient's annual recall appointment, proper fit of the prosthetic restoration or the retention element must be checked. The axial position of the contacts must also be controlled since extra-axial loading may result in loosening of screws and hence in fatigue fractures and fragments may be aspirated or swallowed.

All products labelled ② are intended for single use only. However, if this product is used several times, the following risks may arise: cross contamination, malfunction, inaccurate fit, etc.

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on processing. Any occurring claims for damages may only be made up to the value of our products.

The user and/or patient shall be obliged to report serious incidents that occur in relation to the product to the manufacturer and/or the competent authority of the Member State in which the user and/or patient is registered.

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