crea.lign

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Die layering-technique für removable, fixed and implant-supported restorations.

c) Fabrication of indirect inlays and onlays as well as direct and indirect veneers

d) Customizing or repairing veneers made of composite or ceramic, dentures, crowns and bridges both for direct and indirect indications

3. Contraindication

So far, no contraindications have been reported. If intolerance is suspected, this product may only be used after prior allergological assessment and proof that there is no allergy

4. Hazard and safety information

4.1 crea.lign

Signal word:

Warning

Hazard statements:

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

Precautionary statements:

P280 Wear protective gloves / protective clothing / eye protection.

P302+P350 IF ON SKIN: Gently wash with plenty of soap and water.

P333+P313 If skin irritation or rash occurs: Get medical advice / attention.

P501 Dispose of contents / container in accordance with local / regional / national / international regulations.

Product contains: tetramethylene dimethacrylate

4.2 crea.lign Transpa Clear

Signal word:

Warning

Hazard statements:

H315 Causes skin irritation.

H412 Harmful to aquatic life with long lasting effects.

Precautionary statements:

P280 Wear protective gloves / protective clothing / eye protection.

P302+P350 IF ON SKIN: Gently wash with plenty of soap and water.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P501 Dispose of contents/container in accordance with local / regional / national / international regulations.

4.3 crea.lign Paste

Signal word:

Warning

Hazard statements:

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

Precautionary statements:

P280 Wear protective gloves / protective clothing / eye protection.

P302+P350 IF ON SKIN: Gently wash with plenty of soap and water.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P501 Dispose of contents / container in accordance with local / regional / national / international regulations.

4.4 crea.lign Opaker

Signal word:

Danger

Hazard statements:

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

H412 Harmful to aquatic life with long lasting effects.

Precautionary statements:

P261 Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 Wear protective gloves / protective clothing / eye protection.

P302+P350 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P303+P313 IF ON SKIN: Gently wash with plenty of soap and water.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P501 Dispose of contents / container in accordance with local / regional / national / international regulations.

Product contains: 2-propenoic acid reaction products with pentaerythritol (PETIA/PETTA), diphenyl(2,4,6-trimethylbenzoyl)-phosphine oxide (TPO)

4.5 crea.lign Modelling Liquid

Signal word:

Warning

Hazard statements:

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

Precautionary statements:

P280 Wear protective gloves / protective clothing / eye protection.

P302+P350 IF ON SKIN: Gently wash with plenty of soap and water.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P501 Dispose of contents / container in accordance with local / regional / national / international regulations.

5. Durability and storage

crea.lign (dentine / enamel / incisal / modifier / transpa / gum), crea.lign paste, crea.lign Modelling Liquid and crea.lign opaker must be stored at 4 - 25 ° C and protected against heat and sunlight. These products must not be used after the expiry date.

6. Processing

6.1 Conditioning of metal frameworks

The metal framework is sandblasted with 110 µm aluminium oxide at a pressure of 3 to 4 bars (surface enlargement). Then the surface is wetted with metal adhesive, such as MKZ-Primer (bonding agent for bonding metal / composite); allow the primer to evaporate subsequently (adhere to the processing instructions for the bonding agent). Alternatively, silanizing using inorganic procedures, e.g. Silano Pen, can be carried out. After conditioning the metal, apply a thin coat (max. 0.1 mm) of the respective shade of crea.lign opaker and polymerize (e.g. using bre.Lux Power Unit, Dencatolor XS). To achieve complete coverage of the metal surface, apply at least 2 thin layers of crea.lign opaker and carry out final polymerization for 360 sec. When using mechanical macroretentions, first apply a layer of dual-curing combo.lign opaker and polymerize for 180 sec (adhere to the instructions for use).

Note: After sandblasting, the metal framework should not be cleaned with steam; use a brush (if required, use alcohol) for cleaning.

b) Fabrication of veneer crowns and bridges with metal, ceramic and composite frameworks in the layering technique for removable, fixed and implant-supported restorations.

c) Fabrication of indirect inlays and onlays as well as direct and indirect veneers

d) Customizing or repairing veneers made of composite or ceramic, dentures, crowns and bridges both for direct and indirect indications

3. Contraindication

So far, no contraindications have been reported. If intolerance is suspected, this product may only be used after prior allergological assessment and proof that there is no allergy

4. Hazard and safety information

4.1 crea.lign

Signal word:

Warning

Hazard statements:

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

Precautionary statements:

P280 Wear protective gloves / protective clothing / eye protection.

P302+P350 IF ON SKIN: Gently wash with plenty of soap and water.

P333+P313 If skin irritation or rash occurs: Get medical advice / attention.

P501 Dispose of contents / container in accordance with local / regional / national / international regulations.

Product contains: tetramethylene dimethacrylate

4.2 crea.lign Transpa Clear

Signal word:

Warning

Hazard statements:

H315 Causes skin irritation.

H412 Harmful to aquatic life with long lasting effects.

Precautionary statements:

P280 Wear protective gloves / protective clothing / eye protection.

P302+P350 IF ON SKIN: Gently wash with plenty of soap and water.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P501 Dispose of contents/container in accordance with local / regional / national / international regulations.

4.3 crea.lign Paste

Signal word:

Warning

Hazard statements:

H317 May cause an allergic skin reaction.

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