

Introducing our New Product Line



ARTICULATOR

- 360 degree visibility with posterial incisal guide pin
- Magnetic condyles fixed at 30 degree for free movement articulation



ALU BIG

- 45 minutes denture process using pouring technique
- Clean finish denture using silicone instead of stone
- Adjustable height for any size models
- Three axis locking position



Introductory Offer
10% OFF
All Fast Protec Products
PROMO CODE : LIT-FPINTRO
** Excludes Easy Bar Refills **

Introducing our New Product Line



EASY BAR SYSTEM

- 245 Preformed Bars made of Stainless Steel
- Available in 49 different sizes (from 6 to 30 mm considering center to center distance)
- 1.5 mm thickness
- Bars can be adapted to the different positions of the implants by bending
- Excellent for Implant Impression Jigs
- Measuring Gage
- Easy Bar Template
- Tweezers



** Easy Bar refills available **

Introductory Offer
10% OFF
All Fast Protec Products
PROMO CODE : LIT-FPINTRO
** Excludes Easy Bar Refills **

