

BREDDENT VKS ATTACHMENT RETENTION LEVEL CHART

VKS-SG 1.7mm matrices			
43006680	Level # 1	4 lbs of retained force	Green
43006660	Level # 2	6 lbs of retained force	Yellow
43006640	Level # 3	8 lbs of retained force	Red
430 SG 104 <i>super snap</i>	Level # 4	10 lbs of retained force	Light Green
430 SG 105 <i>super snap</i>	Level # 5	12 lbs of retained force	Light Yellow
430 SG 106 <i>super snap</i>	Level # 6	14 lbs of retained force	Light Red
430 SG 107 <i>super snap</i>	Level # 7	16 lbs of retained force	Very Light Green
430 SG 108 <i>super snap</i>	Level # 8	18 lbs of retained force	Very Light Yellow
430 SG 109 <i>super snap</i>	Level # 9	20 lbs of retained force	Very Light Red
430 SG 110 <i>super snap</i>	Level # 10	22 lbs of retained force	Lightest Green



VKS-SG 2.2mm matrices			
43005410	Level # 1	4 lbs of retained force	Green
43005420	Level # 2	6 lbs of retained force	Yellow
43005430	Level # 3	8 lbs of retained force	Red
430 SG 204 <i>super snap</i>	Level # 4	10 lbs of retained force	Light Green
430 SG 205 <i>super snap</i>	Level # 5	12 lbs of retained force	Light Yellow
430 SG 206 <i>super snap</i>	Level # 6	14 lbs of retained force	Light Red
430 SG 208 <i>super snap</i>	Level # 8	18 lbs of retained force	Very Light Yellow
430 SG 209 <i>super snap</i>	Level # 9	20 lbs of retained force	Very Light Red
430 SG 210 <i>super snap</i>	Level # 10	22 lbs of retained force	Lightest Green

VKS-OC 1.7mm matrices			
43006550	Level # 1	4 lbs of retained force	Green
43006590	Level # 2	6 lbs of retained force	Yellow
43006560	Level # 3	8 lbs of retained force	Red
430 OC 104 <i>super snap</i>	Level # 4	10 lbs of retained force	Light Green
430 OC 105 <i>super snap</i>	Level # 5	12 lbs of retained force	Light Yellow
430 OC 106 <i>super snap</i>	Level # 6	14 lbs of retained force	Light Red
430 OC 107 <i>super snap</i>	Level # 7	16 lbs of retained force	Very Light Green



VKS-OC 2.2mm matrices			
43005440	Level # 1	4 lbs of retained force	Green
43005450	Level # 2	6 lbs of retained force	Yellow
43005460	Level # 3	8 lbs of retained force	Red
430 OC 204 <i>super snap</i>	Level # 4	10 lbs of retained force	Light Green
430 OC 205 <i>super snap</i>	Level # 5	12 lbs of retained force	Light Yellow
430 OC 206 <i>super snap</i>	Level # 6	14 lbs of retained force	Light Red