


BioHPP

High Performance Polymer



Bionics – Surpassed only by nature itself

Catalogue

bredent

Analogue Fabrication

BioHPP granules

Individual dosing, e.g. for crowns, bridges and abutments



Colour	Quantity	REF *	REF
	20 g	540F2PB2	54F2PP01
	100 g	540F2PB3	54F2PP02
 dentin-shade 2	20 g	540F2P22	54F2PP03
	100 g	540F2P23	54F2PP04
 gingiva-shade	20 g	540F2P32	54F2PP05
	100 g	540F2P33	54F2PP06

Not all variants of the products are available in every country!

BioHPP pellets

For major work such as secondary and tertiary treatment.



Colour	Quantity	pce		REF *	REF
	20 g	10	Pellet (15 mm), 2 g	540F2PB6	54F2PP09
	40 g	20	Pellet (15 mm), 2 g	540F2PB7	54F2PP10
	75 g	10	Pellet (25 mm), 7,5 g	540F2PB4	54F2PP07
	150 g	20	Pellet (25 mm), 7,5 g	540F2PB5	54F2PP08
 dentin-shade 2	20 g	10	Pellet (15 mm), 2 g	540F2P26	54F2PP13
	40 g	20	Pellet (15 mm), 2 g	540F2P27	54F2PP14
	75 g	10	Pellet (25 mm), 7,5 g	540F2P24	54F2PP11
	150 g	20	Pellet (25 mm), 7,5 g	540F2P25	54F2PP12
 gingiva-shade	20 g	10	Pellet (15 mm), 2 g		54F2PP17
	40 g	20	Pellet (15 mm), 2 g		54F2PP18
	75 g	10	Pellet (25 mm), 7,5 g		54F2PP15
	150 g	20	Pellet (25 mm), 7,5 g		54F2PP16

* Orders placed from outside the EU



BioHPP dentin-shade 2



BioHPP gingiva-shade and the new smart connect system

Analogue Fabrication

Consumables & accessories for the *for2press* 2 system

Product	Content	REF
Sets		
<i>for2press</i> mould set (muffle system) Ø 16 mm, size 3 4-piece set.	Size 3 mould plate 16 mm, size 3 fixing ring, size 3 metal ring, size 3 fleece liners, 50 units	360F2P16
<i>for2press</i> mould set (muffle system) Ø 26 mm, size 9 4-piece set.	Size 9 mould plate 26 mm, size 9 fixing ring, size 9 metal ring, size 9 fleece liners, 50 units	360F2P26
<i>for2press</i> mould set (muffle system) MP, Ø 26 mm, size 10 9-pce set.	Size 10 mould plate, Ø 26 mm press plunger spacer, Ø 26 mm and 30 mm MP spacer centring aid, height limiter, size 10 metal ring and silicone fixing ring, size 10 fleece liners, 50 units	360F2P27
Mould plate		
<i>for2press</i> mould plate size 3 (16 mm)	Unit	360F2PT1
size 9 (26 mm)	Unit	360F2PT4
MP mould plate size 10 (26 mm)	Unit	360F2PT5
Metal rings		
<i>for2press</i> metal ring for size 3	Unit	360F2PM3
for size 9	Unit	360F2PM9
MP metal ring for size 10	Unit	360F2PM8
Fixing ring for metal ring		
<i>for2press</i> fixing ring for metal ring size 3	Unit	360F2PS3
size 9	Unit	360F2PS9
MP metal ring size 10	Unit	360F2PS8
Fleece liners		
<i>for2press</i> fleece liners 1 mm, Metal ring size 3	50 liners, 143 x 66 x 1 mm	360F2PV3
size 9	50 liners, 227 x 66 x 1 mm	360F2PV9
size 10	50 liners, 264 x 66 x 1 mm	360F2PV1
Disposable press plunger		
<i>for2press</i> disposable press plunger (filler) Ø 16 mm	25 units	570F2P16
(filler) Ø 26 mm	14 units	570F2P26
(filler) Ø 30 mm	14 units	570F2P30
Investment material		
Brevest <i>for2press</i> , also suitable for 3D printing	7.35 kg, 35 x 210 g	570F2PV2
Bresol <i>for2press</i> liquid, standard packaging	1000 ml, standard packaging	520F2PL2
Bresol <i>for2press</i> liquid, frost protection packaging, heat pack	1000 ml, frost protection packaging incl. heat pack and ColdMark	520F2PL3
Set: Brevest <i>for2press</i> , 1 litre Bresol, standard packaging	7.35 kg, 35 x 210 g Brevest <i>for2press</i> 1000 ml Bresol <i>for2press</i> liquid	570F2PV3
Set: Brevest <i>for2press</i> , 1 litre Bresol, frost protection packaging, heat pack	7.35 kg, 35 x 210 g Brevest <i>for2press</i> 1000 ml Bresol <i>for2press</i> liquid Frost protection packaging and heat pack	570F2PV4
Accessories		
perfect2press, channels	30 units	430F2P30
<i>for2press</i> MP spacer, Ø 26 mm, size 10	Unit	360F2PT6
<i>for2press</i> MP centring aid, Ø 26 mm, size 10	Unit	360F2PZ1
<i>for2press</i> MP height limiter, Ø 26 mm, size 10	Unit	360F2PH1
<i>for2press</i> MP spacer Ø 26 mm, size 10	Unit	360F2PA1
<i>for2press</i> MP spacer Ø 30 mm, size 10	Unit	360F2PA2
Material sample BioHPP, 9.5 mm	Unit	992M043D
Digital thermometer	Unit	99300364
Temperature sensor Type K	Unit	99300365
Digital thermometer with temperature sensor	Unit/unit	99300363
Compressed air maintenance unit with connectors	Unit	99300720
Literature		
+ <i>for2press</i> processing instructions	German (available in 5 other languages)	000538GB
+ Operation and maintenance instructions	multilingual	000014EX

Analogue Fabrication

for2press 2

Unit

for2press 2 unit	REF 14000610
+Digital thermometer with temperature sensor	REF 99300363
+Compressed air maintenance unit	REF 99300720
+Operation and maintenance instructions (multilingual)	REF 000014EX
+for2press processing instructions German (available in 5 other languages)	REF 000538GB
+course 50% discount voucher	REF 000537GB

Basic set for beginners

for2press 2 Basic set 8-pce	REF 14000620
enthält	
for2press 2 unit	REF 14000610
BioHPP granulate (20 g)	REF 54F2PP01 REF* 540F2PB2
for2press mould set (muffle system) Ø 16 mm set	REF 360F2P16
for2press filler 16 mm (25 units)	REF 570F2P16
Brevest for2press, 35 x 210 g, 1 l Bresol, standard packaging	REF 570F2PV3
perfect2press, channels (30 units)	REF 430F2P30
BioHPP elaboration and polishing set 15-pce	REF 330F2P10

* Orders placed from outside the EU

Upgrade set

for2press 2 upgrade set 6-pce	REF 14000623
contains	
1. for2press mould set (muffle system) Ø 26 mm, 4 pcs	REF 360F2P26
for2press mould plate size 9, Ø 26 mm	REF 360F2PT4
for2press silicone ring size 9	REF 360F2PR9
for2press metal ring size 9	REF 360F2PM9
for2press fleece liners size 9, 50 pieces	REF 360F2PV9
2. for2press mould Set MP Ø 26 mm + 30 mm, 9 pcs	REF 360F2P27
for2press mould plate size 10	REF 360F2PT5
for2press MP spacer, Ø 26 mm, size 10	REF 360F2PT6
for2press MP spacer, 26 mm	REF 360F2PA1
for2press MP spacer, 30 mm	REF 360F2PA2
for2press MP centring aid	REF 360F2PZ1
for2press MP height limiter	REF 360F2PH1
for2press MP metal ring size 10	REF 360F2PM8
for2press MP silicone fixing ring size 10	REF 360F2PS8
for2press fleece liners size 10, 50 pieces	REF 360F2PV1
3. visio.link PMMA & Composite	REF VLPMMMA10
4. for2press disposable press plunger (filler) Ø 26 mm	REF 570F2P26
5. for2press disposable press plunger (filler) Ø 30 mm	REF 570F2P30
6. BioHPP for2press pellets 75 g (10 pieces, 7,5 g each)	REF 54F2PP07 REF* 540F2PB4

* Orders placed from outside the EU



Not all variants of the products are available in every country!

Digital Fabrication



breCAM.BioHPP

The milling blank

FBioHPP is available as breCAM.BioHPP milling disc for processing in the CAD/CAM workflow. The milling discs can be processed with any open milling system and the

systems from Zirkozahn® and AmannGirrbach® using the specially developed breCAM.cutter.

Product/system	Product/size	REF	Colour
breCAM.BioHPP for open systems	98.5 mm, Ø 12 mm	54002029	white
	98.5 mm, Ø 16 mm	54002030	
	98.5 mm, Ø 20 mm	54002031	
	98.5 mm, Ø 24 mm	54002032	
breCAM.BioHPP for Zirkozahn®	95 mm, Ø 20 mm	54002091	white
breCAM.BioHPP for AmannGirrbach®	84.5 mm, Ø 20 mm	54002111	
breCAM.BioHPP open systems	98.5 mm, Ø 12 mm	54002069	dentine-shade 2
	98.5 mm, Ø 16 mm	54002070	
	98.5 mm, Ø 20 mm	54002071	
	98.5 mm, Ø 24 mm	54002072	
breCAM.BioHPP for Zirkozahn®	95 mm, Ø 20 mm	54002101	dentine-shade 2
	95 mm, Ø 25 mm	54002102	
breCAM.BioHPP for AmannGirrbach®	84.5 mm, Ø 20 mm	54002121	

breCAM.cutter

The patented cutter geometry simplifies dry processing of thermoplastic materials such as BioHPP, PMMA, composites and wax and reduces the reworking time.

for the following machine types	Model	Shaft diameter [mm]	tool diameter [mm]	Number of cutting edges	REF (previous)	REF (new)
AmannGirrbach®	Mikro 4X, Mikro 5X, Mikro IC, Motion 2	3,00	0,60	1	breCAMY28	347B1PN06
AmannGirrbach®	Mikro 4X, Mikro 5X, Mikro IC, Motion 2	3,00	1,00	1	breCAMY32	347B1PN10
AmannGirrbach®	Mikro 4X, Mikro 5X, Mikro IC, Motion 2	3,00	2,50	1	breCAMY31	347B1PN25
Dentsply® - Degudent®	Cercon brain, Cercon brain expert	3,50	1,00	1	breCAMY22	350B1PN10
Dentsply® - Degudent®	Cercon brain, Cercon brain expert	3,50	2,00	1	breCAMY24	350B1PN20
Dentsply® - Degudent®	Cercon brain, Cercon brain expert	3,50	3,00	1	breCAMY26	350B1PN30
Dentsply® - Sirona®	MC X5	3,00	0,50	1	-	342B1PN05
Dentsply® - Sirona®	MC X5	3,00	1,00	1	-	343B1PN10
Dentsply® - Sirona®	MC X5	3,00	2,50	1	-	344B1PN25
Imes-Icore®	Coritec one, Coritec 140i, Coritec 245i, Coritec 250i	3,00	0,60	1	-	348B1PN06
Imes-Icore®	Coritec one, Coritec 140i, Coritec 245i, Coritec 250i	3,00	1,00	1	breCAMX47	348B1PN10
Imes-Icore®	Coritec one, Coritec 140i, Coritec 245i, Coritec 250i	3,00	2,50	1	breCAMX48	348B1PN25
Imes-Icore®	Coritec 350i, Coritec 450i, Coritec 650i, Coritec 750i, Coritec 850i	6,00	0,60	1	breCAMY55	653B1PN06
Imes-Icore®	Coritec 350i, Coritec 450i, Coritec 650i, Coritec 750i, Coritec 850i	6,00	1,00	1	breCAMY57	653B1PN10
Imes-Icore®	Coritec 350i, Coritec 450i, Coritec 650i, Coritec 750i, Coritec 850i	6,00	2,50	1	breCAMY34	653B1PN25
Roland®	DWX-4, DWX-50, DWX-51D, DWX-52D	4,00	1,00	1	breCAMX53	450B1PN10
Roland®	DWX-4, DWX-50, DWX-51D, DWX-52D	4,00	2,00	1	breCAMX54	450B1PN20
Wieland®	Zenotec select/hybrid, Zenotec mini	3,00	1,00	1	breCAMX47	338B1PN10
Wieland®	Zenotec select/hybrid, Zenotec mini	3,00	2,00	1	breCAMX48	338B1PN20
vhf	K3, K4, Z4	3,00	1,00	1	breCAMX67	335B1PN10
vhf	K3, K4, Z4	3,00	2,00	1	breCAMX69	335B1PN20
vhf	S1, S2, K5, R5	3,00	1,00	1	400310KU	340B1PN10
vhf	S1, S2, K5, R5	3,00	2,00	1	400320KU	340B1PN20
Zirkozahn®	M1 Soft Milling Unit	3,00	0,60	1	breCAMY53	357B1PN06
Zirkozahn®	M1 Soft Milling Unit	3,00	1,00	1	breCAMY05	357B1PN10
Zirkozahn®	M1 Soft Milling Unit	3,00	2,00	1	breCAMY07	357B1PN20
Zirkozahn®	M1 Wet, M1 Wet Heavy, M4 Wet Heavy, M5 Wet Heavy, M6 Wet Heavy	6,00	0,50	1	-	650B1PN05
Zirkozahn®	M1 Wet, M1 Wet Heavy, M4 Wet Heavy, M5 Wet Heavy, M6 Wet Heavy	6,00	1,00	1	-	650B1PN10
Zirkozahn®	M1 Wet, M1 Wet Heavy, M4 Wet Heavy, M5 Wet Heavy, M6 Wet Heavy	6,00	2,00	1	-	650B1PN20

Please check the settings of your CAM software before using the cutters.

Indications - analogue


smart connect for BioHPP

The individual screw connection

Thanks to the smart connect for BioHPP screw connection system, the advantages of the physiological BioHPP substructure material are now also available for partially removable restorations. The screw connection can be individually adapted to the restorative situation.

smart connect for BioHPP	REF
smart connect closing ring	540SCB01
smart connect screw short 1.4 mm	540SCB02
smart connect cartridge short M1,4	540SCB03
Fixation screw M 1.4, 2 pieces	36001030
Screw driver short	33000690

smart connect for BioHPP set 5-pce	REF 540F2PSC
contains	
smart connect closing ring	540SCB01
smart connect screw short 1.4 mm	540SCB02
smart connect cartridge short M1,4	540SCB03
Fixation screw M1,4	36001030
Screw driver short	33000690



Indications - analogue

BioHPP titanium base

The conventionally modelled abutment

The main element of the BioHPP elegance hybrid abutments is the titanium base that can be pressed over easily and conveniently with BioHPP using the individual modelling process - the for2press system. Available for the following implant systems:



Company/supplier	System	Product	Titanium grade	REF
bredent medical	SKY classic narrowSKY blueSKY	SKY elegance titanium base incl. SKY screw 2.2	4	SKYETB00
	copaSKY	copaSKY elegance titanium base ML	4	COPAETBM
Straumann®	Bone Level™	elegance titanium base for Straumann® Bone Level™ RC™ 4.1	4	STRFTBR0
		elegance titanium base for Straumann® Bone Level™ NC™ 3.3	4	STRFTBN0
		Screw M 1.6 for Straumann® Bone Level™ RC™ 4.1	Niob	STRSCR16
		Screw M 1.6 for Straumann® Bone Level™ NC™ 3.3	Niob	STRSCN16
	Tissue Level™	elegance titanium base for Straumann® Tissue Level™ RN™ 4.8	5	STTFTRN0
		elegance titanium base for Straumann® Tissue Level™ WN™ 6.5	5	STTFTWN0
Screw for Straumann® Tissue Level™ RN™ 4.8		5	STTSCRNO	
Dentsply®	Astra Tech® OsseoSpeed®	elegance titanium base for Astra Tech® OsseoSpeed® 3.5/4.0	4	AST FTB40
		elegance titanium base for Astra Tech® OsseoSpeed® 4.5/5.0	4	AST FTB50
		Screw M 1.6 for Astra Tech® OsseoSpeed® 3.5/ 4.0	5	ASTSC400
		Screw M 2.0 for Astra Tech® OsseoSpeed® 4.5/ 5.0	5	ASTSC500
	Friadent® XiVe®	elegance titanium base for Dentsply® Friadent® XiVe® 3.4	5	FRIXFTNP
		elegance titanium base for Dentsply® Friadent® XiVe® 3.8	5	FRIXFTRP
		elegance titanium base for Dentsply® Friadent® XiVe® 4.5	5	FRIXFTWP
		Screw for Dentsply® Friadent® XiVe® 3.4	5	FRIXSCNP
		Screw for Dentsply® Friadent® XiVe® 3.8	5	FRIXSCRP
		Screw for Dentsply® Friadent® XiVe® 4.5	5	FRIXSCWP
Nobel Biocare®	Nobel Active™	elegance titanium base for Nobel Active™ RP™ 4.3	4	NBFTBRP0
		elegance titanium base for Nobel Active™ NP™ 3.5	4	NBFTBNP0
		Screw for Nobel Active™ RP™ 4.3	5	NBSCRPO0
		Screw for Nobel Active™ NP™ 3.5	5	NBSCNP00
	Nobel Brånemark™	elegance titanium base for Nobel Brånemark™ NP™ 3.5	5	NMFTBNP0
		elegance titanium base for Nobel Brånemark™ RP™ 4.1	5	NMFTBRP0
		elegance titanium base for Nobel Brånemark™ WP™ 5.1	5	NMFTBWPO
		Screw for Nobel Brånemark™ NP™ 3.5	5	NMSCBNP0
		Screw for Nobel Brånemark™ RP™ 4.3	5	NMSCBRP0
		Screw for Nobel Brånemark™ WP™ 5.0	5	NMSCBWPO
	Nobel Replace™	elegance titanium base for Nobel Replace™ NP™ 3.5	5	NRFTBNP0
		elegance titanium base for Nobel Replace™ RP™ 4.3	5	NRFTBRP0
elegance titanium base for Nobel Replace™ WP™ 5.0		5	NRFTBWPO	
Screw for Nobel Replace™ NP™ 3.5		5	NRSCBNP0	
Camlog®	Screw/Rootline™	elegance titanium base A for Camlog® Screw/Rootline™ 3.8	5	CET3800A
		elegance titanium base B for Camlog® Screw/Rootline™ 3.8	5	CET3800B
		elegance titanium base A for Camlog® Screw/Rootline™ 4.3	5	CET4300A
		elegance titanium base B for Camlog® Screw/Rootline™ 4.3	5	CET4300B
Biohorizons®	Biohorizons®	Screw M 1.6 for Camlog® Screw/Rootline™ 3.8/ 4.3	5	CAMSCM16
		elegance titanium base for Biohorizons® 3.7	5	BIOHFTNP
		elegance titanium base for Biohorizons® 4.25	5	BIOHFTRP
		elegance titanium base for Biohorizons® 5.25	5	BIOHFTWP
		Screw for Biohorizons® 3.7	5	BIOHSCNP
		Screw for Biohorizons® 4.25	5	BIOHSCRP
Sweden & Martina Kohno®	Sweden & Martina Kohno®	Screw for Biohorizons® 5.25	5	BIOHSCWP
		elegance titanium base for Sweden & Martina Kohno® 3.8	5	SWEMFTNP
		elegance titanium base for Sweden & Martina Kohno® 4.25	5	SWEMFTRP
		elegance titanium base for Sweden & Martina Kohno® 5.0	5	SWEMFTWP
		Screw for Sweden & Martina Kohno® 3.8	5	SWEMSCNP
		Screw for Sweden & Martina Kohno® 4.25	5	SWEMSCRP
Zimmer®	Screw Vent®	Screw for Sweden & Martina Kohno® 5.0	5	SWEMSCWP
		elegance titanium base for Zimmer® Screw Vent® 3.5 MIS 3.5/ 4.5	5	ZSVMFTNP
		elegance titanium base for Zimmer® Screw Vent® 4.5	5	ZSVOFTRP
		elegance titanium base for Zimmer® Screw Vent® 5.7 MIS 5.7	5	ZSVMFTWP
		Screw for Zimmer® Screw Vent® 3.5 MIS 3.5/ 4.5	5	ZSVMSCNP
		Screw for Zimmer® Screw Vent® 4.5	5	ZSVOSCRP
3i®	Osseotite®	Screw for Zimmer® Screw Vent® 5.7 MIS 5.7	5	ZSVMSCWP
		elegance titanium base for 3i® Osseotite® 3.4	5	3ICEFTNP
		elegance titanium base for 3i® Osseotite® 4.1	5	3ICEFTRP
		elegance titanium base for 3i® Osseotite® 5.0	5	3ICEFTWP
		Screw for 3i® Osseotite® Certain 3.4/ 4.1/ 5.0	5	3ICESC00

All names marked with ® or ™ are protected brands and/or company names or third-party rights holders.

Indications - analogue

BioHPP elegance hybrid abutment

The prefabricated abutment




The prefabricated abutment after intraoral processing by the dentist.

Photo: Dr. Goldschmid, Lingen, Dentalatelier MDT Martina Brüffer, Osnabrück, Germany

You can find the right shape for most patients from a selection of four different sizes (S, M, L, XL) and two angles (0° and 15°). The prefabricated abutments are designed as a prepared stump, providing the best possible basis for further processing by a doctor or dental technician. If the designs on offer will not work, particularly due to the condition of the gingiva, you can simply order a customised thermopressed abutment from your laboratory.

Customised. Comfortable. Resilient.


No matter if modelled individually, milled or prefabricated - the BioHPP elegance hybrid abutments allow to achieve particularly aesthetic and physiological results.

Company/supplier	System	Product	REF	Colour
bredent medical	SKY classic narrowSKY blueSKY	SKY elegance S Abutment 0°	SKYEES00	 white
		SKY elegance M Abutment 0°	SKYEEM00	
		SKY elegance L Abutment 0°	SKYEEL00	
		SKY elegance S Abutment 15°	SKYEES15	
		SKY elegance M Abutment 15°	SKYEEM15	
		SKY elegance S Abutment 15° R	SKYEES16	
	SKY elegance M Abutment 15° R	SKYEEM16		
	copaSKY	copaSKY elegance M Abutment 0°	COPAEM00	
Straumann®	Bone Level™ RC	elegance S Abutment 0°	STRKSR00	
		elegance M Abutment 0°	STRKMR00	
		elegance L Abutment 0°	STRKLR00	
		Bone Level™ NC	elegance S Abutment 0°	STRKSN00
Dentsply®	Astra OsseoSpeed™	elegance S Abutment 0° 4.5/5.0	ASTKS500	
		elegance M Abutment 0° 4.5/5.0	ASTKM500	
		elegance L Abutment 0° 4.5/5.0	ASTKL500	
		elegance S Abutment 15° 4.5/5.0	ASTKS515	
		elegance M Abutment 15° 4.5/5.0	ASTKM515	
		elegance XL Abutment 0° 4.5/5.0	ASTKXL50	
		elegance S Abutment 0° 3.5/4.0	ASTKS400	
		elegance S Abutment 15° 3.5/4.0	ASTKS415	
Nobel Biocare®	Nobel Active™ RP	elegance S Abutment 0°	NBKSRP00	
		elegance M Abutment 0°	NBKMRP00	
		elegance L Abutment 0°	NBKLRP00	
		elegance XL Abutment 0°	NBKXLRP0	
		elegance S Abutment 15°	NBKSRP15	
		elegance M Abutment 15°	NBKMRP15	
		Nobel Active™ NP	elegance S Abutment 0°	NBKSNP00
			elegance S Abutment 15°	NBKSNP15

Indications - digital

BioHPP elegance prefab

The abutment for digital fabrication in the CAD/CAM workflow for SKY classic / blueSKY / narrowSKY and other standard implant systems.



Company/ supplier	System	Product	Ø Implant connection	Titanium grade	REF	REF scan abutment	Colour
bredent medical	SKY classic, blueSKY, narrowSKY	BioHPP SKY elegance prefab Set	3.5 mm	4	SKYEPFST	SKYUSCAI	 dentine-shade 2
3i Biomet®	Osseotite Certain®	BioHPP elegance prefab	3.4 mm	5	BOC34B201	OC34SCAN	
			4.1 mm	5	BOC41B201	OC41SCAN	
			5.0 mm	5	BOC50B201	OC50SCAN	
Dentsply®	Astra OsseoSpeed®	BioHPP elegance prefab	3.5 mm	5	BOS35B201	OS35SCAN	
			4.5 mm	5	BOS45B201	OS45SCAN	
	Friadent® XiVE®	BioHPP elegance prefab	3.4 mm	5	BXV34B201	XV34SCAN	
			3.8 mm	5	BXV38B201	XV38SCAN	
			4.5 mm	5	BXV45B201	XV45SCAN	
Biohorizons®	Biohorizons Internal®	BioHPP elegance prefab	3.7 mm	5	BBI37B201	BI37SCAN	
			4.25 mm	5	BBI42B201	BI42SCAN	
			5.25 mm	5	BBI52B201	BI52SCAN	
Camlog®	CONELOG®	BioHPP elegance prefab	3.3 mm	4	BCL33B201	CL33SCAN	
			3.8 mm	4	BCL38B201	CL38SCAN	
			4.3 mm	4	BCL43B201	CL43SCAN	
			5.0 mm	4	BCL50B201	CL50SCAN	
	SCREW/ROOT- LINE	BioHPP elegance prefab	3.8 mm	5	BSL38B201	SL38SCAN	
			4.3 mm	5	BSL43B201	SL43SCAN	
Nobel Biocare®	Brånemark TiUnite® MkIII	BioHPP elegance prefab	3.5 mm	5	BTU35B201	TU35SCAN	
			4.1 mm	5	BTU41B201	TU41SCAN	
			5.1 mm	5	BTU51B201	TU51SCAN	
	NobelActive	BioHPP elegance prefab	3.5 mm	5	BNA35B201	NA35SCAN	
			4.3 mm	5	BNA43B201	NA43SCAN	
	NobelReplace® Tapered/ NobelParallel™ CC/ NobelReplace® CC	BioHPP elegance prefab	3.5 mm	5	BRP35B201	RP35SCAN	
			4.3 mm	5	BRP43B201	RP43SCAN	
			5.0 mm	5	BRP50B201	RP50SCAN	
	Straumann®	Bone Level™	BioHPP elegance prefab	3.3 mm	4	BBL33B201	BL33SCAN
4.1 mm				4	BBL41B201	BL41SCAN	
synOcta®		BioHPP elegance prefab	4.8 mm	4	BSO48B201	SO48SCAN	
			6.5 mm	4	BSO65B201	SO65SCAN	
Sweden & Martina®	Kohno	BioHPP elegance prefab	3.8 mm	5	BKH38B201	KH38SCAN	
			4.25 mm	5	BKH42B201	KH42SCAN	
			5.0 mm	5	BKH50B201	KH50SCAN	
Zimmer®	Tapered Screw Vent®	BioHPP elegance prefab	3.5 mm	5	BSV35B201	SV35SCAN	
			4.5 mm	5	BSV45B201	SV45SCAN	
			5.7 mm	5	BSV57B201	SV57SCAN	

All names marked with ® or ™ are protected brands and/or company names or third-party rights holders.

Accessories

Finishing and polishing

Avoid plaque accumulation with the polishing protocol for the BioHPP finishing and polishing set.



BioHPP finishing and polishing set, 15 pcs.		REF 330F2P10
contains		
Diamond disc Giflex-TR Master x-tray		34000M25
TC cutter 1.4 mm - grenade		H272M814
TC cutter 4.0 mm - grenade		H274M840
TC cutter 6.0 mm - grenade		H274M860
TC cutter 2.3 mm - pear		H237M823
TC cutter 2.3 mm - conical round		H200M823
TC cutter 4.0 mm - inverted conical		H263M840
Spiral mandrel		00303023
Screwable mandrel		34000660
Round brushes goat-hair white, double the bristles, Ø 22 mm		35000550
Ceragum coarse, Cylinder		PWKG0612
Ceragum coarse, Lens		PLKG2212
Cotton buff with shaft, Ø 22 mm		35000650
Acrypol polishing paste, 170 g		52000170
Abraso-Starglanz high luster polishing paste, 50 g		52000163



Use the BioHPP telescope set to achieve optimum friction values for your telescopic restorations.



BioHPP telescope set, 8 pcs.		REF 330F2P11
contains		
TC cutter 1.0 mm - cylinder, round		F137M810
TC cutter 1.5 mm - cylinder, round		F137M815
TC cutter 2.3 mm - cylinder, round		F137M823
Pen-shaped brushes Chungking black, 2 mm		35000430
Pen-shaped brushes Chungking black, 4 mm		35100410
Pen-shaped brushes goat-hair white, 2 mm		35000440
Pen-shaped brushes goat-hair white, 4 mm		35000420
Acrypol polishing paste, 170 g		52000170



Photo: MDT Jens-Christian Fehsenfeld, DENTEC Dentallabor, Münster, Germany



visio.link


The perfect combination of composite and resins. visio.link is a universal, light-curing PMMA- and composite primer and ensures an optimum bond between PMMA-based denture resins, high-impact PMMA composite materials (novo.lign veneers/neo.lign denture teeth/acrylic teeth) and combo.lign.



novo.lign

Veneers

The novo.lign veneers feature an extended neck and a natural depth thanks to the morphological layer structure. The gold standard in implant prosthetics.



Quick Reference Card


The practical guide for orientation when using BioHPP:

- Preparation and impression taking
- Sterilising, conditioning and bonding
- Cleaning in the practice and laboratory as well as chairside

bre dent group portfolio




Other offers that may be of interest to you:



REF 000534GB



REF 000547GB



REF 000538GB

000535GB-2021.013 Mistake and subject to change reserved

