

LOCATOR® Attachment Components

Implant Abutments Compatibility Chart

LOCATOR SYSTEM COMPONENTS				PRODUCT CODE			PRICE (USD)					
LOCATOR Implant Abutment				(see table below)			\$169.99					
LOCATOR Core Tool (includes abutment driver)				8393			\$104.99					
LOCATOR Male Processing Package (2 Pack)				8519-2			\$99.99					
LOCATOR Male Processing Package (10 Pack)				8519-10			\$449.99					
LOCATOR Extended Range Male Processing Package (2 Pack)				8540-2			\$99.99					
LOCATOR Extended Range Male Processing Package (10 Pack)				8540-10			\$449.99					
LOCATOR IMPLANT ABUTMENT PRODUCT CODES												
IMPLANT DESCRIPTION			HEIGHT OF ABUTMENT TISSUE CUFF ABOVE THE IMPLANT									
			0mm	0.73mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm	7.0mm	8.0mm
a IMPLANT												
4.0				1421 (.5mm)		1422	1423	1424	1425			
AB DENTAL (888.823.5553)												
3.75 NP			8661		8662	8663 (2.5mm)	8664 (3.5mm)	8665 (4.5mm)	8625 (5.5mm)	8626 (6.5mm)		
4.50 WP			8671		8672	8673	8674	8675	8627	8628		
ACE SURGICAL (800.441.3100)												
4.0, 5.0 (4.1) EXTERNAL HEX				8681		8682	8683	8684	8685	8696		
4.0 & 5.0 INTERNAL CONNECTION				ACE SURGICAL (Available from ACE SURGICAL)								
ALMITECH (516.343.3050)												
3.3, 4.1, 4.8			1438		1439	1440	1441	1442	1443			
ALPHA BIO TEC LTD												
3.75				1192 (0.3mm)	1193	1194	1195	1196	1197	1198		
AMERROSS DENTAL IMPLANT SYSTEM / FORMALLY DICOA (1.877.342.6211)												
LOCATOR Implant Abutments (Available from AMERROSS)												
ANTHOGRYR												
3.4 SHORT						1697						
3.4 MEDIUM									1698			
3.4 TALL										1699		
ASTRA TECH / DENTSPLY IH AB (800.531.3481)												
3.5, ST (M2.0 thread)				8636 (0.81mm)		8637	8638	8639	8640			
4.0, ST (M2.0 thread)				8656 (0.5mm)		8657	8658	8659	8660			
4.5, 5.0, ST (M2.0 thread) (Lilac)				8686 (0.5mm)		8687	8688	8689	8690			
3.5/4.0 (M1.6 thread) (Aqua)					8115	8116	8117	8118	8119			
BEGO IMPLANT SYSTEMS												
LOCATOR Implant Abutments available from BEGO.												
BICON (800.882.4266)												
2.0 POST					8746	8747	8748	8749	8750			
3.0 POST					8930	8931	8932	8933	8934			
BIOHORIZONS (888.246.8338)												
3.5 (yellow) EXTERNAL HEX				8531 (1.25mm)	8532	8533	8534	8535				
4.0 (green) EXTERNAL HEX				8541	8542	8543	8544	8545				
5.0 (blue) EXTERNAL HEX				8551 (1.5mm)	8552 (2.5mm)	8553 (3.5mm)	8554 (4.5mm)	8555 (5.5mm)				
3.0 INTERNAL LASER LOCK				1927	1928	1929	1930					
3.5 (yellow) INTERNAL HEX			8661	8662	8663 (2.5mm)	8664 (3.5mm)	8665 (4.5mm)	8625 (5.5mm)	8626 (6.5mm)			
4.5 (green) INTERNAL HEX			8671	8672	8673	8674	8675	8627	8628			
5.7 (blue) INTERNAL HEX				8190	8191	8192	8193	8194	8195			
BIO-LOK (800.789.0830)												
3.45 (green) INTERNAL AND EXTERNAL HEX					8561	8562	8563	8564	8565			
4.0 (yellow) INTERNAL AND EXTERNAL HEX				8681		8682	8683	8684	8685	8696		
5.0 (red) INTERNAL AND EXTERNAL HEX					8581 (1.4mm)	8582	8583	8584	8585			
BIOMET (3i) (800.342.5454)												
3.25 EXTERNAL HEX					8441	8442	8443	8444	8445	8446		
3.75 MINIPLANT EXTERNAL HEX				8681		8682	8683	8684	8685	8686		
4.0 CYLINDRICAL EXTERNAL HEX				8681		8682	8683	8684	8685	8686		
5.0 EXTERNAL HEX				8820 (0.6mm)		8821	8822	8823	8824			
TG					8621	8622	8623	8624				
3.4 CERTAIN INTERNAL MICRO IMPLANT					8431	8432	8433	8434	8435	8436		
4.1 CERTAIN INTERNAL (Fits Certain 5.0 & 6.0)				8630		8631	8632	8633	8634	8635		
BIOTEC SRL												
3.25 TRE						8441	8442	8443	8444	8445		8446
3.75 TRE				8681		8682	8683	8684	8685	8686		8687
3.75 SEI			8090			8091	8092	8093	8094	8095		8096
8 NECK						8950	8951	8952	8953	8954		8955
BIOTECH												
KONTACT												
BLUE SKY Bio (718.376.0422)												
3.7 INTERNAL HEX			8661		8662	8663 (2.5mm)	8664 (3.5mm)	8665 (4.5mm)	8625 (5.5mm)	8626 (6.5mm)		
4.7 INTERNAL HEX			8671		8672	8673	8674	8675	8627	8628		
3.5mm INTERNAL TRILOBE			8756		8757	8758	8759	8760	8734	8735		
4.3mm INTERNAL TRILOBE			8761		8762	8763	8764	8765	8753	8754		
5.0mm INTERNAL TRILOBE					8766	8767	8768	8769	8770			
ONE-STAGE: REGULAR INTERNAL				8950		8951	8952	8953	8954	8955		
BREDENT												
4.1 Sky					1827	1828	1829	1830	1831	1832		
BTI - BIOTECHNOLOGY INSTITUTE SL												
4.1 EXTERNAL HEX					8681	8682	8683	8684	8685	8686		
4.1 INTERNAL					8291	8292	8293	8294	8295	8296		
CAMLOG (877.537.8862)												
LOCATOR® Implant Abutments (Available from Camlog Biotechnologies, Henry Schien USA/Canada)												
COWELL MEDI USA, Inc. (714.52.1800)												
ATLAS Sub.1					1735	1736	1737	1738	1739	1740		
DENTAURUM												
3.4, 3.7 and 4.5 TIOLOX												
LOCATOR Implant Abutments (Available from DENTAURUM)												
DENTEGRIS / DENTAL TECH												
3.75			8052		8053	8054	8055	8056	8057			

LOCATOR® Attachment Components

Implant Abutments Compatibility Chart

IMPLANT DESCRIPTION	HEIGHT OF ABUTMENT TISSUE CUFF ABOVE THE IMPLANT									
	0mm	0.73mm 8058 (0.63 mm)	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm	7.0mm	8.0mm
4.5				8059	8060	8061	8062	8063		
DENTIUM USA (877.304.6752)										
3.4/3.8, 4.3/4.8, 4.8W IMPLANTUM		8656 (0.5mm)		8657	8658	8659	8660			
DENTSPLY IMPLANTS MFG. (DENTSPLY FRIADENT) (800.531.3481)										
FRIALIT-2 & XIVE 3.4			8805	8806	8807	8808	8809			
FRIALIT-2 & XIVE 3.8			8810	8811	8812	8813	8814			
FRIALIT-2 & XIVE 4.5			8815	8816	8817	8818	8819			
FRIALIT 2- & XIVE 5.5			8139 (1.3mm)	8140	8141	8142	8143			
ANKYLOS C/X (Fits Ankylos Plus)				1729	1730	1731	1732	1733		
DSI										
3.8 SM (NARROW) DIO			1351	1352	1353	1354	1355	1356		
4.5/5.3 SM (REGULAR) DIO			1361	1362	1363	1364	1365	1366		
Dyna Dental										
4.0			8999							
Equinox										
6.0		1192 (0.3mm)		1193	1194	1195	1196	1197		
FAIR IMPLANT GmbH										
4.2 Fair			1948 (1.5mm)	1949	1950	1951	1952	1953		
FMZ GmbH ALPHATECH										
LOCATOR Implant Abutments (Available from FMZ-GmbH Alphatech)										
GC JAPAN										
LOCATOR Implant Abutments (Available from GC Japan)										
GEASS										
LOCATOR Implant Abutments (Available from Geass)										
GLOBAL IMPLANT SOLUTIONS, LLC (781.232.9477)										
4.0 OSSEOLINK		1910 (0.2mm)	1911	1912	1913	1914	1915			
IMPLADENT (PHIBO)										
TSH S2			1916	1917	1918	1919	1920	1921		
TSH S3-S4			1942	1943	1944	1945	1946	1947		
TSH S5				1922	1923	1924	1925	1926		
TSA S3			1755	1756	1757	1758	1759	1760		
TSA S4			1761 (1.4mm)	1762	1763	1764	1765	1766		
TSA S5				1767	1768	1769	1770	1771		
IMTEC CORPORATION (800.879.9799)										
SELF TAPPING 3.75 EXTERNAL HEX		8681		8682	8683	8684	8685	8686		
UNIVERSAL 3.4, 4.0 EXTERNAL HEX		8681		8682	8683	8684	8685	8686		
MINI 3.3 EXTERNAL HEX		8856 (0.91mm)		8857	8858	8859	8860			
3.5, 4.3 ENDURE		8994 (0.63mm)		8995	8996	8997	8998			
IMZ										
IMZ 3.3 EXTERNAL HEX		8681		8682	8683	8684	8685	8686		
IMZ 3.3 NON-HEX FLAT TOP	8020		8021	8022	8023	8024	8025			
IMZ 4.0 EXTERNAL HEX		8681		8682	8683	8684	8685	8686		
IMZ 4.0 NON-HEX FLAT TOP	8030 (0.35mm)		8031	8032	8033	8034	8035			
INTRA-LOCK SYSTEM INTERNATIONAL (877.330.0338)										
3.75, 4.0 THREADED EXTERNAL HEX		8681		8682	8683	8684	8685	8686		
4.3 THREADED CONIC EXTERNAL HEX			8781	8782	8783	8784	8785			
4.75 THREADED EXTERNAL HEX			8727 (1.42mm)	8728	8729	8730	8731			
5.0 THREADED CONIC EXTERNAL HEX			8727 (1.42mm)	8728	8729	8730	8731			
INTERNAL CONNECT LOCATOR® Implant Abutments (Available from Intra-Lock Corporation)										
KAT IMPLANTS, LLC (603-427-0084)										
LOCATOR® Implant Abutments (Available from Kat Implants)										
KEYSTONE/LIFECORE BIOMEDICAL (866.902.9272)										
3.8		LOCATOR Implant Abutments available from KEYSTONE DENTAL INC.								
4.8										
6.5										
RENOVA 3.75	8090		8091	8092	8093	8094	8095	8096		
RENOVA 4.5, 4.75		8125 (0.63mm)		8126	8127	8128	8129	8130		
RESTORE 3.3 and 3.4 (SD) EXTERNAL HEX			8196	8197	8198	8199	8200			
RESTORE 3.75, 4.0 (RD) EXTERNAL HEX		8681		8682	8683	8684	8685	8686		
RESTORE 5.0, 6.0 (WD) EXTERNAL HEX			8727 (1.42mm)	8728	8729	8730	8731			
SUSTAIN 4.2 EXTERNAL HEX		8681		8682	8683	8684	8685	8686		
SUSTAIN 4.0, 4.7	8651		8652	8653	8654	8655				
STAGE-1		8950		8951	8952	8953	8954	8955		
PRIMA CONNEX LOCATOR Implant Abutments (Available from Keystone Dental)										
KLOCKNER										
LOCATOR Implant Abutments (Available from SOADCO, S.L.)										
KYOCERA										
4.2 POI		8716 (0.35mm)	8717	8718	8719	8720	8721			
LASAK										
3.7	1207		1208	1209	1210	1211	1212			
MEDENTIKA										
4.0, 4.5, 5.5 M-Implant			1624	1625	1626	1627	1628	1629		
MEDENTIS MEDICAL GmbH										
2.0 ICX			1643	1644	1645	1646	1647	1648		
MIS MEDICAL IMPLANT SYSTEM (866.797.1333)										
3.75, 4.2 INTERNAL (SP) IMPLANT	8661		8662	8663 (2.5mm)	8664 (3.5mm)	8665 (4.5mm)	8625 (5.5mm)	8626 (6.5mm)		
5.0 INTERNAL (WP) IMPLANT	8671		8672	8673	8674	8675	8627	8628		
EXTERNAL HEX (SP) IMPLANT		8681		8682	8683	8684	8685	8686		
EXTERNAL HEX (WP) IMPLANT		8820 (0.6mm)		8821	8822	8823	8824			
M & K DENTAL GmbH										
LOCATOR Implant Abutments available from M & K GmbH										

LOCATOR® Attachment Components

Implant Abutments Compatibility Chart

IMPLANT DESCRIPTION	HEIGHT OF ABUTMENT TISSUE CUFF ABOVE THE IMPLANT									
	0mm	0.73mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm	7.0mm	8.0mm
NEODENT CONE MORSE CONNECTION										
LOCATOR ATTACHMENT SYSTEM										
CM PLATFORM				2168	2169	2170	2171	2172		
LOCATOR R-TX ATTACHMENT SYSTEM										
CM PLATFORM			32301-01	32301-02	32301-03	32301-04	32301-05	32301-06		
LOCATOR F-TX ATTACHMENT SYSTEM										
CM PLATFORM			12301-02	12301-03	12301-04	12301-05	12301-06			
NEOSS LIMITED										
4.0			8451	8452	8453	8454	8455	8456		
NOBEL BIOACRE (800.993.8100)										
NOBELACTIVE (NP) 3.5 Conical Connection			2069	2070	2071	2072	2073	2074		
NOBELACTIVE (RP) 4.3/5.0 Conical Connection			2075	2076	2077	2078	2079	2080		
SELECT, REPLACE 3.5 (Purple) Internal Tri Lobe	8756		8757	8758	8759	8760	8734	8735		
SELECT, REPLACE 4.3 (Gold) Internal Tri Lobe	8761		8762	8763	8764	8765	8753	8754		
SELECT, REPLACE 5.0 (Blue) Internal Tri Lobe			8766	8767	8768	8769	8770			
NOBEL SPEEDY REPLACE (GROOVY) 4.0	8761		8762	8763	8764	8765	8753	8754		
BRANEMARK 3.3 (NP) EXTERNAL HEX		8105		8106	8107	8108	8109			
BRANEMARK 3.75, 4.0, 5.0 (RP) EXTERNAL HEX		8681		8682	8683	8684	8685	8696		
BRANEMARK 5.0 (WP) EXTERNAL HEX			8825 (1.1mm)	8826	8827	8828	8829			
NOBEL SPEEDY EXT HEX (GROOVY) 4.0		8681		8682	8683	8684	8685	8696		
NOBEL SPEEDY EXT HEX (GROOVY) 3.5		8105		8106	8107	8108	8109			
REPLACE 3.5 EXTERNAL HEX			8676	8677	8678	8679	8680			
REPLACE 4.3 EXTERNAL HEX			8781	8782	8783	8784	8785			
REPLACE 5.0 EXTERNAL HEX			8835 (1.37mm)	8836	8837	8838	8839			
HEX-LOCK (HL) 3.25 EXTERNAL HEX			8646 (1.5mm)	8647	8648	8649	8650			
HEX-LOCK (HL) 3.8, 4.5 EXTERNAL HEX		8681		8682	8683	8684	8685	8696		
HEX-LOCK (HL) 5.0 EXTERNAL HEX			8835 (1.37mm)	8836	8837	8838	8839			
STERIOSS 3.25 NON-HEXED			(There is no LOCATOR Implant Abutment available for this product)							
STERIOSS 3.8 NON-HEXED	8651		8652	8653	8654	8655				
OCO BIOMEDICAL, INC. (800.228.0477)										
3.25 TSI		1892 (.70mm)		1893	1894	1895	1896		1897	
4.0 TSI		8950		8951	8953	8954	8955			
5.0 TSI				1898	1899	1900	1901		1902	
OSSTEM / HIOSSEN USA										
4.8 SS INTERNAL/HG REGULAR		8950		8951	8952	8953	8954	8955		
3.5 GS MINI INTERNAL/HG MINI			1561	1562	1563	1564	1565	1566		
4.0 GS STANDARD INTERNAL/HG STANDARD			1551	1552	1553	1554	1555	1556		
3.3/3.5 US MINI EXTERNAL HEX				8106	8107	8108	8109			
3.75/4.0 US REGULAR/HU REGULAR EXT HEX		8681		8682	8683	8684	8685	8696		
OSTEO-IMPLANT CORPORATION (800.654.5560)										
3.25 SD, 3.75 ST, 4.0 DT		8681		8682	8683	8684	8685	8696		
PERIOSEAL INC (888.446.7126)										
LOCATOR® Implant Abutments (Available from Perioseal Inc.)										
PRODENT ITALIA										
LOCATOR Implant Abutments (Available from PRODENT ITALIA)										
PROWITAL GmbH (WIELAND)										
LOCATOR Implant Abutments (Available from Prowital)										
RIEMSER ARZNEIMITTEL (CURASAN) (888.237.2767)										
REVOIS 3.8, 4.3 AND 5.0			1213	1214	1215	1216	1217	1218		
REVOIS II		1303 (0.2mm)	1304	1305	1306	1307	1308			
SARGON (888.688.6684)										
3.8 EXTERNAL HEX			8722 (1.8mm)		8723	8724	8725	8726		
SCHUTZ DENTAL										
IMPLA LOCATOR Implant Abutments (Available from SCHUTZ DENTAL GROUP)										
SERF IMPLANTER INNOVATION										
LOCATOR Implant Abutments (Available from SERF Implanter Innovation)										
SIC INVENT AG										
LOCATOR Implant Abutments (Available from SIC Invent GmbH)										
SIMPACT USA (973.588.8930)										
LOCATOR Implant Abutments (Available from Simpac)										
SOUTHERN IMPLANT USA (866.700.2100)										
3.5 INTERNAL TRI-LOBE	8756		8757	8758	8759	8760	8734	8735		
4.3 INTERNAL TRI-LOBE	8761		8762	8763	8764	8765	8753	8754		
5.0/6.0 INTERNAL TRI-LOBE			8766	8767	8768	8769	8770			
4.8 OCT (ITST) INTERNAL OCTAGON		8950		8951	8952	8953	8954	8955		
STERNGOLD IMPLAMED (800.243.9942)										
HEX CYLINDER 3.3, 4.0 EXTERNAL HEX		8681		8682	8683	8684	8685	8696		
HEX SCREW 3.75, 4.0, 5.0 EXTERNAL HEX		8681		8682	8683	8684	8685	8696		
STRAUMANN (800.448.8168)										
LOCATOR ATTACHMENT SYSTEM										
STRAUMANN TISSUE LEVEL CONNECTION - NNC			2056	2057	2058	2059	2060	2061		
STRAUMANN TISSUE LEVEL CONNECTION - RN			8742	8743	8744	8745	8751	8752		
STRAUMANN TISSUE LEVEL CONNECTION - WN			8870	8871	8872	8873	8874			
STRAUMANN BONE LEVEL CONNECTION - NC				8422	8423	8424	8425	8426		
STRAUMANN BONE LEVEL CONNECTION - RC			8170	8171	8172	8173	8174	8175		
LOCATOR R-TX ATTACHMENT SYSTEM										
STRAUMANN TISSUE LEVEL CONNECTION - NNC			30405-01	30405-02	30405-03	30405-04	30405-05	30405-06		
STRAUMANN TISSUE LEVEL CONNECTION - RN			30406-01	30406-02	30406-03	30406-04	30406-05	30406-06		
STRAUMANN TISSUE LEVEL CONNECTION - WN			30407-01 (1.5mm)	30407-02	30407-03	30407-04	30407-05	30407-06		
STRAUMANN BONE LEVEL CONNECTION - NC			30401-01	30401-02	30401-03	30401-04	30401-05	30401-06		
STRAUMANN BONE LEVEL CONNECTION - RC			30402-01	30402-02	30402-03	30402-04	30402-05	30402-06		
LOCATOR F-TX ATTACHMENT SYSTEM										
STRAUMANN TISSUE LEVEL CONNECTION - NNC			10405-02	10405-03	10405-04	10405-05	10405-06			
STRAUMANN TISSUE LEVEL CONNECTION - RN			10406-02	10406-03	10406-04	10406-05	10406-06			
STRAUMANN TISSUE LEVEL CONNECTION - WN			10407-02 (2.5mm)	10407-03	10407-04	10407-05	10407-06			
STRAUMANN BONE LEVEL CONNECTION - NC			10401-01	10401-02	10401-03	10401-04	10401-05	10401-06		
STRAUMANN BONE LEVEL CONNECTION - RC			10402-01	10402-02	10402-03	10402-04	10402-05	10402-06		
LOCATOR Implant Abutments (Available from Straumann)										
SUDIMPLANT										
Sudimplant / TBR	LOCATOR Implant Abutments (Available from Sudimplant)									

LOCATOR® Attachment Components

Implant Abutments Compatibility Chart

IMPLANT DESCRIPTION	HEIGHT OF ABUTMENT TISSUE CUFF ABOVE THE IMPLANT									
	0mm	0.73mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm	7.0mm	8.0mm
SWEDEN & MARTINA										
3.8 PREMIUM KOHNO			1675	1676	1677	1678	1679	1680		
4.25 PREMIUM KOHNO			1681	1682	1683	1684	1685	1686		
SYBRON / IMPLANT DIRECT (www.implantdirect.com)										
3.7 LEGACY (3.5 Collar)	8661		8662	8663 (2.5mm)	8664 (3.5mm)	8665 (4.5mm)	8625 (5.5mm)	8626 (6.5mm)		
4.7 LEGACY (4.5 Collar)	8671		8672	8673	8674	8675	8627	8628		
5.7 LEGACY (5.5 Collar)			1471	1472	1473	1474	1475	1476		
3.5 RePLUS	8756		8757	8758	8759	8760	8734	8735		
4.3 RePLUS	8761		8762	8763	8764	8765	8753	8754		
5.0 RePLUS			8766	8767	8768	8769	8770			
3.7 SCREWPLANT			1371	1372	1373	1374	1375	1376		
4.7 SCREWPLANT			1381 (1.5mm)	1382	1383	1384	1385	1386		
5.7 SCREWPLANT				1391	1392	1393	1394	1395		
3.7 SID				1250	1251	1252	1253	1254	1255	
3.7 SWISSPLANT	1821 (0.4mm)		1822	1823	1824	1825	1826			
SYBRON / SYBRON IMPLANT SOLUTIONS (800.537.7824)										
4.1 ENDOPORE (.7mm EXTERNAL HEX)		8681		8682	8683	8684	8685	8696		
3.25 ENTEGRA (.7mm EXTERNAL HEX)		8681		8682	8683	8684	8685	8696		
4.0 ENTEGRA (.7mm EXTERNAL HEX)		8681		8682	8683	8684	8685	8696		
3.5 ENDOPORE (.9mm EXTERNAL HEX)			8851	8852	8853	8854	8855			
3.25 ENTEGRA (.9mm EXTERNAL HEX)			8851	8852	8853	8854	8855			
4.1 ENDOPORE (.9mm EXTERNAL HEX)			(There is no LOCATOR Implant Abutment available for this product)							
4.7 ENTEGRA (.9mm EXTERNAL HEX)			8691 (1.35mm)	8692	8693	8694	8695			
5.0 ENDOPORE (.9mm EXTERNAL HEX)			8691 (1.35mm)	8692	8693	8694	8695			
4.1, 5.0 ENDOPORE ONE-STAGE		8950		8951	8952	8953	8954	8955		
4.1 ENDOPORE (INTERNAL HEX)	8145		8146	8147	8148	8149	8150			
5.0 ENDOPORE (INTERNAL HEX)			8164	8165	8166	8167	8168	8169		
3.25 / 3.3 PITT-EASY			8311	8312	8313	8314	8315	8316		
3.75 / 3.8 PITT-EASY			8321	8322	8323	8324	8325	8326		
4.0 / 4.1 PITT-EASY			8331	8332	8333	8334	8335	8336		
4.9 PITT-EASY			8341	8342	8343	8344	8345	8346		
3.5 PRO SERIES (XRT BONE LEVEL)			1631	1632	1633	1634	1635	1636		
4.1 PRO SERIES (XRT BONE LEVEL)			1637	1638	1639	1640	1641	1642		
4.8 PRO SERIES (XRT BONE LEVEL)		8950		8951	8952	8953	8954	8955		
4.1 / 4.8 PRO SERIES (TL TISSUE LEVEL)		8950		8951	8952	8953	8954	8955		
3.5 SYBRON EXTERNAL HEX PLUS		8105		8106	8107	8108	8109			
4.1 SYBRON EXTERNAL HEX PLUS		8681		8682	8683	8684	8685	8696		
5.1 SYBRON EXTERNAL HEX PLUS (WP)			8825 (1.1mm)	8826	8827	8828	8829			
TATUM SURGICAL (888.360.5550)										
LOCATOR® Implant Abutments (Available from Tatum Surgical)										
TEKKA										
3.2 KRESTAL			1691	1692	1693	1694	1695	1696		
3.4 TEKKA	1220		1221	1222	1223	1224	1225			
4.0 TEKKA			1226	1227	1228	1229	1230	1231		
4.5 TEKKA			1232	1233	1234	1235	1236	1237		
5.0 TEKKA			1238	1239	1240	1241	1242	1243		
6.0 TEKKA			1244	1245	1246	1247	1248	1249		
I-KONE		1937 (0.7mm)	1938 (1.5mm)		1939		1940		1941	
THOMMEN MEDICAL (866.319.9800)										
SPI LOCATOR Implant Abutments (Available from Thommen Medical)										
ZIMMER - CALCITEK (800.854.6691)										
INTEGRAL 3.25 SD	8611		8612	8613	8614	8615				
INTEGRAL 4.0	8651		8652	8653	8654	8655				
OMNIOLOC 3.25	8611		8612	8613	8614	8615				
OMNIOLOC 4.0	8666		8667	8668	8669	8670				
SPLINE 3.25			8701	8702	8703	8704	8705			
SPLINE 3.75, 4.0			8711	8712	8713	8714	8715	8697	8698	8699
SPLINE 5.0		8736 (0.86mm)		8737	8738	8739	8740			
SWEDE-VENT 4.0 EXTERNAL HEX		8681		8682	8683	8684	8685	8696		
SWEDE-VENT TL EXTERNAL HEX		8681		8682	8683	8684	8685	8696		
TAPERLOCK EXTERNAL HEX		8681		8682	8683	8684	8685	8696		
THREADLOCK 4.0 EXTERNAL HEX		8681		8682	8683	8684	8685	8696		
ZIMMER - PARAGON (800.854.6691)										
3.5mm DIAMETER COLLAR: (green)										
SCREW-VENT 3.3, 3.7, 4.1	8661		8662	8663 (2.5mm)	8664 (3.5mm)	8665 (4.5mm)	8625 (5.5mm)	8626 (6.5mm)		
4.5mm DIAMETER COLLAR: (purple)										
SCREW-VENT 4.7	8671		8672	8673	8674	8675	8627	8628		
5.5mm DIAMETER COLLAR: (gold)										
SCREW-VENT 5.7			8180	8181	8182	8183	8184	8185		
SWISSPLUS 3.8		8935 (0.75mm)		8936	8937	8938				
SWISSPLUS 4.8			8939 (1.29mm)	8940	8941	8942				
ADVENT 4.5		8841 (0.71mm)		8842	8843	8844	8845	8846		
ZITERION										
Locator Implant Abutments (Available from Ziterion)										
ZL - MICRODENT										
ZL DURAPLANT LOCATOR Implant Abutments (Available from ZL-Microdent GmbH & Co. KG)										